Proteinase K

Code No. 9033
Size : 5 ml

Shipping at – 20°C
Store at – 20°C

Lot No.

Concentration : 20 mg/ml

Volume : 100 mg

Expiration Date :

Description :

Proteinase K is a serine protease with high activity. This enzyme activity is stable under the existence of calcium salt and increases when supplied with protein denaturant such as SDS or urea.

This enzyme can be available for a wide range of substrate. Especially it degrades esters and peptide bonds preferentially next to C-termini of hydrophobic, sulfuric or aromatic amino acid.

Storage Buffer :

20 mM Tris-HCl, pH7.4
1 mM CaCl₂
50% Glycerol

Source : *Tritirachium album*

Unit definition :

One unit of enzyme activity is defined as the amount required to liberate Folin positive amino acid corresponded to 1 μ mol of tyrosine using bovine hemoglobin denatured by urea, at 37°C for 1 minute.

Specific Activity :

> 30 U/mg protein

Purity :

Nuclease activity was not detected in any of the following tests, as judged from the agarose gel electrophoresis pattern :

- 1) After incubation of 1 μ g of λ DNA-*Hind* III fragments with 1 μ l of this enzyme for 24 hours at 37°C.
- 2) After incubation of 1 μ g of supercoiled pBR322 DNA with 1 μ l of this enzyme for 1 hour at 37°C.
- 3) After incubation of 1 μ g of 16S or 23S rRNA with 1 μ l of this enzyme for 24 hours at 37°C.

Applications :

- Purification of DNA, RNA and phage
- Preparation of chromosomal DNA for *in situ* hybridization, fingerprinting, colony hybridization, plaque hybridization and pulsed-field electrophoresis

Application example :

0.01 M Tris-HCl, pH7.8
0.01 M EDTA
0.5% SDS

50 - 500 μ g/ml Proteinase K

↓

37 - 56°C, 1hr - overnight

↓

Use the mixture for the next reaction, following phenol/chloroform extraction and ethanol precipitation.

References :

- 1) Ebeling, W., *et al.* (1974) *Eur. J. Biochem.*, **47**, 91-97.
- 2) Sambrook, J., *et al.* (1989) *Molecular Cloning: A Laboratory Manual*, 2nd edition, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- 3) Gunkel, F. A. and Gassen, H. G. (1989) *Eur. J. Biochem.* **179**, 185-194.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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