EASY Dilution (for Real Time PCR)

Code No. 9160

Size : 1 ml × 8 tubes

Shipping at - 20°C

Stored at - 20°C

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Lot No.

Description:

This product is used to dilute template DNA or RNA when preparing the serial dilution of standard for establishing the standard curve in real time PCR. When template DNA and RNA are diluted with sterilized water or TE buffer, they can stick to the inside tube which may prevent correct dilution. This product helps to avoid occurring this absorption and the correct dilution at lower concentrations is possible. In addition, as this product does not include nucleic acid like tRNA or rRNA, non-specific amplification derived from those RNA will not happen.

Note:

EASY Dilution is designed to use with real time PCR products produced by TaKaRa Bio Inc. The combination of this product and other company's products hasn't been tested.

Protocol:

The concentration of this product is × 1. It can be directly used for dilution of template DNA or RNA. Serial dilution of the standard solution is added into reverse transcription or PCR reaction at an arbitrary amount. But it is recommended to add a certain volume in a series of experiments.

(Example)

Preparation for serial dilution of cDNA* for drawing the standard curve

Prepare the serial dilution series of cDNA at the concentration of 5 ng/μl, 500 pg/μl, 50 pg/μl, 5 pg/μl and 0.5 pg/μl by diluting a portion of the cDNA solution (50 ng/μl) by 10-fold each.

*: The cDNA was synthesized from total RNA through reverse transcription.

By adding 2 μl of serial dilution to each reaction, the standard curve toward the cDNA in the range of 1 pg to 100 ng is drawn. (For negative control reaction, add this product by 2 μl each.)

Note:

It is recommended to store this product in dispensing into aliquots convenient for each customer in order to avoid contamination among samples.

Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at www.takara-bio.com.

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