TaKaRa Ex Taq[®] Hot Start Version

Code No. RR006A
Size : 250 units

Shipping at – 20°C
Store at – 20°C

Components :

TaKaRa Ex Taq HS	250 U
10X Ex Taq Buffer (20 mM Mg ²⁺ plus)	1 ml
dNTP Mixture (2.5 mM each)	800 µl

Conc. : 5 units/µl
Volume : 50 µl

Description :

TaKaRa Ex Taq HS is designed for hot start PCR; it includes a neutralizing monoclonal antibody that recognizes Taq DNA polymerase. This antibody inhibits polymerase activity by binding to Taq, thereby preventing nonspecific amplification due to mispriming and/or formation of primer dimers during reaction set-up before thermal cycling. Antibody-mediated repression is released during the initial DNA denaturation step of PCR. This enzyme can be used with standard PCR conditions.

Storage Buffer :

20 mM	Tris-HCl (pH8.0)
100 mM	KCl
0.1 mM	EDTA
1 mM	DTT
0.5%	Tween20
0.5%	Nonidet P-40
50%	Glycerol

Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM	TAPS (pH9.3 at 25°C)
50 mM	KCl
2 mM	MgCl ₂
0.1 mM	DTT
200 µM ea.	dATP·dGTP·dCTP
100 µM	[³ H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

Purity :

Nicking, endonuclease, and exonuclease activity were not detected after incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ DNA, or 0.6 µg of λ-Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

PCR products :

As most PCR products amplified with TaKaRa Ex Taq HS have one A at the 3'-termini, the obtained PCR products can be directly cloned into a T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Application : For DNA amplification by hot start PCR

PCR test :

Good performance in PCR was confirmed by amplification of a 17.5 kb fragment of the β-globin gene from human genomic DNA, and amplification of a 20 kb fragment from λ DNA.

General reaction mixture for PCR (50 µl reaction volume)

TaKaRa Ex Taq HS (5 units/µl)	0.25 µl
10X Ex Taq Buffer	5 µl
dNTP Mixture (2.5 mM each)	4 µl
Template	< 500 ng
Primer 1	10 - 50 pmol (final conc. 0.2 - 1.0 µM)
Primer 2	10 - 50 pmol (final conc. 0.2 - 1.0 µM)
Sterile distilled water	up to 50 µl

(Note) Reaction mixtures can be prepared at room temperature. Keep all reagents on ice.

PCR conditions :

This enzyme can be used with standard PCR conditions, since the monoclonal antibody is denatured during the initial DNA-denaturation step. There is no need for an additional step to denature the anti-Ta_q antibody.

Example : Amplification of 1 kb DNA fragment

98°C	10 sec.] 30 cycles	98°C	10 sec.] 30 cycles
55°C	30 sec.		68°C	1 min.	
72°C	1 min.				

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec. at 98°C, or 20 - 30 sec. at 94°C is recommended.

TaKaRa Ex Taq is a registered trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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