mRNA Interferase™ -MazF

Code No. 2415A **Size:** **1,000 U**
Conc.: **20 U/μl**

Supplied Reagent:
5X MazF Buffer **1 ml**

Description:

MazF is a toxin protein in the toxin-antitoxin module of *E. coli*, having endoribonuclease activity which specifically cleaves single-stranded RNA at 5' end of ACA sequence. This enzyme does not cleave double-stranded RNA, double-stranded DNA or single-stranded DNA. This product, mRNA Interferase-MazF, is supplied as a fusion protein of *E. coli* MazF and trigger factor which is one of *E. coli* chaperon proteins.

Storage Buffer:

20 mM Sodium phosphate, pH 6.0
0.01% Tween 20
50% Glycerol

Storage: -20°C

Source: Expressed as a recombinant protein in *Escherichia coli*

Unit definition:

One unit is the amount of enzyme that cleaves 1 pmol of standard substrate (ROX-5'-GATAUACATATCT-eclipse; the underlined region are composed of RNA) for 10 min at 37°C and pH 7.5.

Buffer composition for unit definition:

40 mM Sodium phosphate, pH 7.5
0.01% Tween 20

Supplied reagent: 5X MazF Buffer

200 mM Sodium phosphate, pH 7.5
0.05% Tween 20

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Application: Limited digestion of ssRNA

Application Example:

1) Cleavage of 30 base ssRNA including one ACA sequence

30 base ssRNA* (100 pmol/μl)	1.0 μl
mRNA Interferase-MazF (20 U/μl)	1.0 μl
5X MazF Buffer	2.0 μl
DEPC-treated water	6.0 μl
Total	10.0 μl

* RNA sequence: 5'-UAAGAAGGAGAUUAACAAUGAAUCAAUC-3'

The ssRNA was completely cleaved into 16 bases and 14 bases after incubation at 37°C for 10 min.

2) Cleavage of 584 base ssRNA including ACA sequence at 13 sites

584 base ssRNA (1 pmol/μl)	0.5 μl
mRNA Interferase-MazF (20 U/μl)	1.0 μl
5X MazF Buffer	1.0 μl
RNase-free water	2.5 μl
Total	5.0 μl

The ssRNA was cleaved into each fragment after the incubation for 15 min at 37°C.

Note:

1. This enzyme specifically cleaves at ACA sequence of ssRNA, but it may cleave AC site depending on the sequences around the AC site.
2. This enzyme does not cleave dsRNA. Accordingly it may not cleave at all ACA sequences included in RNA due to the influence of the secondary structure of the substrate.
3. The addition of salt, especially divalent cation, e.g. Mg²⁺, into the reaction mixture inhibits the enzyme activity. DTT does not have any effect on the reaction.
4. The cleaved fragment by this enzyme cannot be applied directly to ligate by RNA ligase, because it has 2', 3'-cyclic phosphate at one of its termini and 5'-OH at the other.

References:

- 1) Zhang Y, Zhang J, Hara H, Kato I, and Inouye M. *J Biol Chem.* (2005) **280**: 3143-3150.
- 2) Zhang Y, Zhang J, Hoeflich K P, Ikura M, Qing G, and Inouye M. *Mol Cell.* (2005) **12**: 913-923.
- 3) Munoz-Gomez A J, Santos-Sierra S, Berzal-Herranz A, Lemonnier M, and Diaz-Orejas R. *FEBS Lett.* (2004) **567**: 316-320.

Interferase is a trademark of Takara Bio Inc.

1. PURCHASER shall not transfer or sell copies of PRODUCT, component of PRODUCT or derivatives of PRODUCT, or products obtained through the use of PRODUCT to any other party.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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