pHY300PLK DNA

Code No. 3060 **Size:** **10 µg**
Conc.: 300 - 800 µg/ml

* 2 years from date of receipt under proper storage conditions.

Supplied Reagent:
Host strain *Bacillus subtilis* ISW 1214 (glycerol stock)

100 µl

Host strain is shipped separately.

Host strain *Bacillus subtilis* ISW 1214

Form : glycerol stock (50% glycerol)

Storage : -80°C

Genotype : *hsr M1, leuA8, metB5, tet^s*

Description:

pHY300PLK is a shuttle vector functioning in *E. coli* and *B. subtilis*, which is constructed with a plasmid pACYC177 of *E. coli* and DNA derived from plasmid pAM α 1 of *Streptococcus faecalis*.

This DNA can transform both *E. coli* and *B. subtilis*. pHY300PLK has ampicillin resistant gene and tetracycline resistant gene as selective markers. This vector expresses both genes in *E. coli*, but only tetracycline resistant gene in *B. subtilis*.

Form:

10 mM Tris-HCl, pH 7.4
0.1 mM EDTA

Storage: -20°C

Preparation:

Purification of cccDNA by CsCl-EtBr ultracentrifugation

Base pairs: 4,870 bp

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.


Transformation Efficiency:

E. coli ≥ 2 × 10⁵ / µg DNA

B. subtilis ≥ 2 × 10⁷ / µg DNA*

* The transformation efficiency by electroporation.

Vector map of pHY300PLK:

References:

- 1) Anagnostopoulos C and Spizizen J. *J Bacteriol.* (1961) **81**: 741-746.
- 2) Ikawa S, Shibata T, Matsumoto K, Iijima T, Saito H, and Ando T. *Mol Gen Genet.* (1981) **183**: 1-6.
- 3) Sadaie Y and Kada T. *J Bacteriol.* (1983) **153**: 813-821.
- 4) Ishiwa H and Tsuchida N. *Gene.* (1984) **32**: 129-134.
- 5) Ishiwa H and Shibahara H. *Jpn J Genet.* (1985) **60**: 235-243.

Manufactured by Yakult Co., LTD

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

pHY300PLK DNA

Code No. 3060

<Transformation procedure into *Bacillus subtilis* ISW1214>

In case of the transformation into *B. subtilis* by Method A using competent cells, the plasmid DNA should be amplified with *E. coli recA*⁺ cells, as the transformation efficiency is much dependent on the ratio of dimer plasmid.

On the other hand, in case of the transformation into *B. subtilis* by electroporation (Method B), monomer plasmid amplified with *E. coli recA*⁻ cells is effective for transformation.

[Method A] Transformation using competent cells

(1) Preparation of Competent cell

SP I medium		SP II medium	
(NH ₄)SO ₄	0.2%	(NH ₄)SO ₄	0.2%
K ₂ HPO ₄	1.4%	K ₂ HPO ₄	1.4%
KH ₂ PO ₄	0.6%	KH ₂ PO ₄	0.6%
Na-Citrate-2H ₂ O	0.1%	Na-Citrate-2H ₂ O	0.1%
MgSO ₄ ·7H ₂ O	0.02%	MgSO ₄ ·7H ₂ O	0.02%
*Glucose	0.5%	*Glucose	0.5%
*Casamino acid	0.02%	*Yeast extract	0.02%
*Yeast extract	0.1%	*MgCl ₂	5 mM
*L-Leucine	50 μg/ml	*L-Leucine	5 μg/ml
*L-Methionine	50 μg/ml	*L-Methionine	5 μg/ml

*They must be autoclave individually.

1. First Culture

- Inoculate *B. subtilis* ISW 1214 in 5 ml of L-broth and culture at 37°C with shaking overnight.
- The following morning, inoculate the overnight culture of ISW1214 into fresh SP1 medium (inoculation volume is about 1%) and culture at 37°C with shaking.
- When cells have grown to the end of log phase, add glycerol to the culture to the final concentration of 12.5%.
- Dispense a 400 μl of aliquot of the culture into microcentrifuge tubes and freeze quickly in a ethanol-dry ice bath.
- Keep the frozen cells at -70°C until use.

2. Second culture

- Liquefy the frozen culture at 37°C.
- Dilute it to 7.5X volume with SP1I medium.
- Culture it for 90 minutes at 37°C with shaking. Use this culture of ISW1214 as competent cells.

(2) Transformation

- Mix 50 μl of competent cells solution and 50 μl of plasmids solution (10 - 100 ng DNA).
- Shake at 37°C for 30 minutes.
- Add 100 μl of L-broth, and shake for 60 minutes at 37°C.
- Spread the aliquot on L-plate containing 20 μg/ml tetracycline.*1
- Incubate the plate at 37°C overnight.

[Method B] Transformation by electroporation

(1) Preparation of cells for electroporation

- Inoculate *B. subtilis* ISW1214 in L-broth and culture at 37°C with shaking overnight.
- The following morning, inoculate 2 ml of preculture into 32 ml of fresh medium (L-broth+0.5 M sorbitol) and culture at 37°C with shaking.
- When cells have grown to A₆₆₀=0.85 - 0.95, stop the culture. (about 2.5 hours)
- Keep the culture on ice for 10 min.
- Centrifuge at 5,000 *g* for 5 min.
- Wash the cells 4 times with chilled Solution A (0.5 M sorbitol, 0.5 M mannitol, 10% glycerol).
- Suspend the cells with 0.8 ml of Solution A.
- Dispense in 60 μl each and store at -80°C. Use these cells for electroporation.

(2) Electroporation

- Thaw the cells (60 μl) on ice and add 1 μl of DNA solution (50 ng/ml).
- Transfer the mixture of cells and DNA into a chilled cuvette (0.1 cm) and keep it for 1 - 1.5 min on ice.
- Pulse*2 it, add 1 ml of Solution B (L-broth + 0.5 M sorbitol + 0.38 M mannitol) and culture at 37°C for 3 hours.
- Spread the aliquot on L-plate containing 30 μg/ml tetracycline*1 and incubate the plate at 37°C overnight.

Notes:

- *1 The ampicillin resistant gene of pHY300PLK can be expressed in *E. coli*, but cannot be expressed in *B. subtilis*. So the antibiotics resistance gene markers is the tetracycline resistant gene only in *B. subtilis*.
- *2 The pulse condition of Gene Pulser (BIO-RAD Inc.): 2 kV, 25 μF, 200 Ω (Time constant=4.5 - 5.0)

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Host strain *Bacillus subtilis* ISW 1214 (glycerol stock)

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