# pHSG398 DNA

Code No. 3398

Size: 25  $\mu$ g

Conc.: 0.5  $\mu$ g/ $\mu$ l

\* 2 years from date of receipt under proper storage conditions.

## Description:

pHSG398 is a pUC type of cloning vector plasmid possessing chloramphenicol resistancy. This plasmid is suitable for DNA sequencing by dideoxy method. It may allow successful cloning of DNA fragment longer than the case of M13 phage vector used. Since it has multi-cloning sites on *lacZ* region, any DNA insertion into this plasmid can be easily verified using plates containing IPTG and X-Gal. Moreover, expression of the cloned DNA can be performed by using of *lac* promoter. M13 primers are available for DNA sequencing.

Using Deletion Kit for Kilo-Sequencing (Cat. #6030), DNA sequencing of the cloned DNA with kilo-bases can also be performed.

pHSG398 contains the consensus multi-cloning sites as same as pUC18.

**Form:** 10 mM Tris-HCl, pH 8.0  
1 mM EDTA

**Storage:** -20°C

**Preparation:** cccDNA was purified by CsCl-EtBr ultracentrifugation.

**Base pairs:** 2,227 bp

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.



## Usage:

- Cloning of a target gene and its expression using *lac* promoter
- DNA sequencing by using M13 primers
- Determination of long DNA sequences using Deletion Kit for Kilo-Sequencing (Cat. #6030)

## Reference:

Takeshita S, Sato M, Toba M, Masahash W, and Hashimoto-Gotoh T. *Gene*. (1987) **61**:63-74.

## Multicloning site of pHSG398:

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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