Protein Molecular Weight Marker (High)

Code No. 3451

Size: for 200 lanes

Conc.: 10 $\mu\text{g}/\mu\text{l}$

Supplied Reagents:

5X Loading Buffer

1 ml

1 M DTT (Dithiothreitol)

100 μl

Description:

This marker consists of 5 kinds of proteins with a Molecular Weight range between 44.3 kDa and 200 kDa. Each protein is proportioned to yield uniform band intensity when stained with Coomassie Brilliant Blue R-250 after run on SDS-polyacrylamide gel.

5 μl of 20-fold diluted marker is loaded per lane of SDS-PAGE minigel for standard use. In this case, this marker is sufficient for 200 lanes.

Form:

50 mM	Tris-HCl, pH 6.8
1 mM	EDTA
200 mM	NaCl
50 %	Glycerol

Component Proteins:

Protein	Source	MW (Da)
Myosin	Pig	200,000
β -galactosidase	<i>E. coli</i>	116,000
Phosphorylase B	Rabbit muscle	97,200
Serum Albumin	Bovine	66,409
Ovalbumin	Hen egg White	44,287

Storage:

Protein MW Marker (High), 1 M DTT :	-20°C
5X Loading Buffer :	Store at RT after used.

5X Loading Buffer (Store at RT after used):

200 mM	Tris-HCl, pH 6.8
10%	SDS
0.05%	BPB
50%	Glycerol

Application Example:

- Combine the following reagents in a tube.


1 M DTT	2 μl
5X Loading Buffer	20 μl
- Prepare 20-fold diluted marker by adding the following components to the solution prepared in 1.

Protein MW Marker (High)	5 μl
Sterile purified water	73 μl

* 20-fold diluted marker is stable for 2-3 months for -20°C. It is recommended to store in aliquots for several uses, not to repeat freeze-thaw cycles.
- Mix 20-fold diluted marker well, and heat at 100°C for 5 minutes. Load 5 μl per lane of SDS-PAGE minigel. Run 7.5 - 10.5% SDS-polyacrylamide gel electrophoresis.

- Perform staining with Coomassie Brilliant Blue R-250.

Electrophoresis Result (7.5% SDS-PAGE)

Note:

It is recommended to use 7.5 - 10.5% SDS-polyacrylamide gel for electrophoresis. If the concentration of gel is too high, high molecular weight bands may not separate well or they may not migrate.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Protein Molecular Weight Marker (High)

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