# Glycopeptidase F

(Peptide: *N*-glycosidase F)

Code No. 4450

Size: 25 mU/ vial

## Supplied Reagent :

(A) Denature Buffer	500 $\mu$ l
(B) Native Buffer	500 $\mu$ l
(C) Stabilizer Solution	500 $\mu$ l
(D) Control glycoprotein	10 $\mu$ l

## Common name :

Peptide-*N*<sup>4</sup>- (*N*-acetyl- $\beta$ -glucosaminy) asparagine amidase

Enzyme code: 3.5.1.52

## Source :

*Escherichia coli* carrying the plasmid containing Glycopeptidase F gene

Conc.: 500 mU/ ml

## Form :

Solution in 20 mM sodium phosphate (pH 7.2), 50 mM EDTA and 0.05% sodium azide.

## Substrate specificity :

Specific cleavage of *N*-glycans (GluNAc-Asn bonds) in glycopeptides or glycoproteins. In the case Asn-linked GlcNAc contains  $\alpha$ -1,3 fucosidic linkage, this enzyme do not cleavage the *N*-glycan.

## Definition of activity :

One unit is the amount of enzyme required to hydrolyze 1  $\mu$  mol of dansyl fetuin glycopeptide within 1 min at 37°C, pH 8.5.

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Properties :

- Molecular weight : approx. 36,500 (SDS-PAGE)

### (Note)

Depending on proteolytic activity of host bacteria, the recombinant Glycopeptidase F (GPF) may have 15 or 18 additional amino acid residues at its amino terminal (1). However, its activity on glycoproteins was confirmed to be identical to that of native GPF, by the experiments using RNase B and Fetuin as its substrates.

- Optimum pH : pH 8.6

Storage: -20°C

Once thawed, this product should be stored at 4 °C.  
Avoid freeze-thaw cycles.

## Protocol :

There are 2 ways of De-*N*-glycosylation of Glycoprotein by GPF, native condition and denature condition. Please choose the either native or denature condition that is suited to the purpose of each experiment.

## Supplied Reagents :

(A) Denature Buffer	: 1% SDS/1 M Tris-HCl (pH 8.6)
(B) Native Buffer	: 1 M Tris-HCl (pH 8.6)
(C) Stabilizer Solution	: 5% NP-40
(D) Control glycoprotein	: 10 mg/ml Bovine fetuin

## Sample :

The buffer concentration of samples should be lower than 50 mM. If it is higher than 50 mM, change the buffer concentration less than 50 mM by appropriate method (ex. dialysis) or adjust pH to nearly 8.6 (optimum pH of GPF) by addition of acid or alkaline. It is recommended to check pH when adding sample into buffer to make sure.

Samples should not contain citrate buffer and urea.

## Reaction Condition :

GPF is useful for 3 purposes ; **1.** to check whether the protein sample has *N*-Glycan or not, **2.** to examine the effect of removing *N*-Glycan from protein sample for its activity, and **3.** to remove *N*-Glycan from protein sample and analyze the structure of sugar chain. The reaction condition should be chosen suited for each purpose.

### For purpose 1:

Use 1 mU of GPF per 25  $\mu$ g of glycoprotein and incubate overnight under denature condition.

### For purpose 2:

Incubate under native condition. Under native condition, the cleavage efficiency depends on the kind of glycoprotein. At first, try to use 10 mU of GPF and incubate overnight. If not cleaved, add GPF further and incubate overnight.

### For purpose 3:

Under denature condition, sugar chain is completely removed and it costs the small amount of GPF, but it is necessary to remove detergent. If avoid removing detergent, perform the reaction under native condition while checking the cleavage.

When the amount of glycoprotein is small, the amount of the enzyme may be decreased, but it is not recommended to decrease the amount of GPF less than 0.5 mU per 1 reaction (25  $\mu$ l).

## SDS - PAGE :

For Coomassie staining, the cleavage can be detected when 1 - 5  $\mu$ g of the sample is loaded per lane. For silver staining, the cleavage can be detected when 10 - 50 ng of the sample is loaded per lane.

Choose the appropriate concentration of acrylamide gel for the molecular weight of glycoprotein sample. 10% acrylamide gel is appropriate for control glycoprotein, fetuin. If the content of *N*-Glycan of glycoprotein is less than 5%, it may be necessary to extend the electrophoresis time, as it is difficult to detect the change of molecular weight by the cleavage of sugar chain.

## Note :

This product contains sodium azide as preservative. Sodium azide reacts with many heavy metals to form explosive compounds. It is critical that a large amount of water be used when disposing.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# Glycopeptidase F

(Peptide: *N*-glycosidase F)

Code No. 4450

容量： 25 mU/ vial

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(A) Denature Buffer	500 $\mu$ l
(B) Native Buffer	500 $\mu$ l
(C) Stabilizer Solution	500 $\mu$ l
(D) Control glycoprotein	10 $\mu$ l

● 酵素名 Peptide- $N^4$ - (*N*-acetyl- $\beta$ -glucosaminy) asparagine amidase

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*Escherichia coli* carrying the plasmid containing Glycopeptidase F gene

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## Procedure :

1. Digest Control glycoprotein (D) and samples according to the following flow chart. (Use purified water instead of GPF and perform the same reaction as Negative control.)
2. Analyze the reaction mixture by SDS-PAGE. Fetuin [Control glycoprotein] binds to 3 *N*-Glycan per 1 molecule, so when it is deglycosylated, the band is shifted by SDS-PAGE. (Since fetuin binds to 3 or 4 *O*-glycan in addition to 3 *N*-Glycan, the apparent molecular weight becomes higher than accurate molecular weight.

### 1) Denature Condition

- Prepare 2.5  $\mu$ l of sample (glycoprotein, about 25  $\mu$ g) or (D)
- 2.5  $\mu$ l of Denature Buffer (A) with 0.2 M 2-mercaptoethanol
  - heat at 100°C for 3 min
  - add 5  $\mu$ l of Stabilizer Solution (C)
  - mix
  - add 13  $\mu$ l of purified water
  - mix
  - add 2  $\mu$ l (1 mU) of Glycopeptidase F
  - mix
  - incubate at 37°C for 15 - 20 hours
- Analyze for SDS-PAGE (10% acrylamide gel)

### 2) Native Condition

- Prepare 2.5  $\mu$ l of sample (glycoprotein, about 25  $\mu$ g) or (D)
- 2.5  $\mu$ l of Native Buffer (B)
  - add 20  $\mu$ l (10 mU) of Glycopeptidase F
  - mix
  - incubate at 37°C for 15 - 40 hours
- Analyze for SDS-PAGE (10% acrylamide gel)

## Application : (Experiment with the native enzyme)

Glycoprotein (human  $\alpha$ 1-AGP, bovine fetuin, hen ovalbumin, human transferrin, human IgG ; Table 1) and human plasma were cleaved by GPF both under denature and native condition, and analyzed the removal of sugar chain by SDS-PAGE (10% acrylamide gel).

Table 1

Glycoprotein	Source	molecular weight	Number of <i>N</i> -Glycan
$\alpha$ 1-AGP	human	44,100	5
fetuin	bovine	48,400	3
ovalbumin	hen	45,000	1
transferrin	human	75,000	2
IgG	human	150,000 (H chain, 50,000)	2 (1/mol H chain)

## Result :

Under denature condition, when 25  $\mu$ g of glycoprotein was incubated with 1 mU of GPF for 17 hours, *N*-Glycan was removed from all of glycoprotein completely (Figure 1, lane 2, 6, 10, 14, 17). On the other hand, under native condition, there was a remarkable difference in the cleavage efficiency between the kinds of glycoprotein.

Transferrin and IgG were cleaved completely with 10 mU of GPF for 17 hours (Figure 1, lane 15, 18), but  $\alpha$ 1-AGP, fetuin and ovalbumin were not cleaved completely (Figure 1, lane 3, 7, 11). So it is proved that the extension of incubation time more than 17 hours had no effect on the cleavage.

Human plasma was cleaved completely under both denature and native condition (Figure 1, lane 20, 21, 22). It was observed that some bands except the band corresponded to transferrin and IgG were shifted.

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