# $\alpha$ 2,3-Sialidase

**Code No. 4455**                      **Size :**            **5 U/ vial**  
**Conc. :**                              **50 U/ml**

**Note :**

The storage temperature of this product changed from -20°C to 4°C (Lot. N201AA-).

**Systematic name :** Acylneuraminyl hydrolase

**Enzyme code :** 3. 2. 1. 18

**Storage :**                            4°C

**Source :**

Cloned from *Salmonella typhimurium* LT2 and expressed in *Escherichia coli*.

**Form :**

Solution in 50 mM potassium phosphate (pH6.8), 0.1 M sodium chloride, and 0.02% sodium azide.

**Feature :**

This product is the cloned enzyme. It does not require divalent cation (Ca<sup>2+</sup>, Mg<sup>2+</sup>) or stabilizer (BSA etc.)<sup>1,2</sup>.

**Properties :**

Molecular weight : 41,300

Michaelis constant : Km=0.25 mM (4MUNeu5Ac)  
Km=2.66 mM (PA<sup>3</sup>- 3'-Sialyllactose)

Inhibitors : 2-deoxy-2,3-didehydro-N-acetylneuraminic acid  
Ki=0.38 mM

Stabilizer: Not necessary

Optimum pH : pH 5.5 - 7.0 (Tris-maleate buffer)  
pH 5.5 - 6.0 (potassium phosphate buffer, sodium acetate buffer, Tris-HCl buffer)

**Definition of activity :**

One unit is the amount of enzyme required to hydrolyze 1  $\mu$  mol of 3'-sialyllactose within 1 minute at 37°C at pH 5.5.

**Substrate specificity :**

- This sialidase has about 200-fold kinetic preference of the  $\alpha$  2, 3 sialyl linkage to the  $\alpha$  2, 6 sialyl linkage.<sup>1)</sup>
- The enzyme efficiently cleaves sialyl groups from glycoproteins and glycolipids.
- The enzyme can hydrolyze glycolipids without detergent.
- The enzyme may not release sialic acid from O-linked oligosaccharides on glycoproteins

**Quality Control Data :**

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

**References :**

- 1) Hoyer L L, Roggentin P, Schauer R, and Vimr E R. *J Biochem.* (1991) **110**: 462.
- 2) Hoyer L L, Hamilton A C, Steenbergen S M, and Vimr E R. *Mol Microbiol.* (1992) **6**: 873.
- 3) Warren L. *J Biol Chem.* (1959) **234**: 1971.
- 4) Kondo A, Suzuki J, Kuraya N, Hase S, Kato I, and Ikenaka T. *Agric Biol Chem.* (1990) **54**: 2169.

**Note**

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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阻害剤： 2-deoxy-2,3-didehydro-N-acetylneuraminic acid

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



The reaction mixture is as follows:

<Reaction mixture>

2.1 mM 4MUNeu5Ac	] 3 $\mu$ l
50 mM Sodium acetate buffer (pH 5.5)	
100 mM NaCl	
Enzyme solution	3 $\mu$ l
Total volume	6 $\mu$ l

The reaction mixture is incubated for 10 minutes at 37°C, then 1.5 ml of 0.1 M glycine-NaOH buffer, pH 10.3 is added to 2.5  $\mu$ l of the reaction mixture to stop the reaction.

The amount of 4-MU released during the reaction is measured by using fluorescence photometer (Ex: 350 nm; Em: 448 nm).

## Method 2:

10  $\mu$ l of the enzyme solution is added to 90  $\mu$ l of the substrate in 0.1 M sodium acetate at pH 5.5. The reaction mixture is incubated for 10 minutes at 37°C.

The amount of released sialic acid is measured by TBA method<sup>3)</sup>.

## Method 3:

The reaction mixture is as follows:

<Reaction mixture>

10 $\mu$ M PA sugar chain <sup>4)</sup>	2 $\mu$ l
0.1 M Sodium acetate buffer (pH 5.5)	5 $\mu$ l
Enzyme solution	3 $\mu$ l
Total volume	10 $\mu$ l

The reaction mixture is incubated for 10 minutes at 37°C, then analyzed by HPLC.

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