Lacto-*N*-biosidase from *Streptomyces* sp. 142

Code No. 4456

Size: 100 μ U

Note:

The storage temperature of enzyme has been changed to -20°C from 4°C. (Lot. N106AB -)

Description:

This product specifically hydrolyzes oligosaccharides with type I chain and produces lacto-*N*-biose (Gal β 1-3GlcNAc) but does not hydrolyze oligosaccharides with type II chain. So it can be used to judge the structure of type I, II glycoprotein and glycolipid sugar chain. Also as it can distinguish Lewis^s structure from Lewis^x structure in conjunction use with α -1,3/4-L-Fucosidase (Cat. #4453), this enzyme is useful in analyzing structure and function of glycoprotein and glycolipid sugar-chain. One vial of this product (100 μ l) allows 50 - 100 reactions of enzyme digestion of sugar chain up to 10 pmol.

Source: *Streptomyces* sp.142

Form: Solution in 50 mM sodium acetate buffer, pH 5.5, containing 0.05% Brij-58.

Storage: -20°C

Concentration: 1 μ U/ μ l

Reaction:

Specific hydrolysis of oligosaccharides with type I chain (Gal β 1-3GlcNAc) and production of lacto-*N*-biose (Gal β 1-3GlcNAc).

Definition of activity:

One unit is the amount of enzyme required to hydrolyze 1 μ mol of pyridylamino lacto-*N*-tetraose within 1 minute at 37°C at pH 5.5.

Assay of activity:

Incubated with 2 μ M PA-lacto-*N*-tetraose (PA-Sugar Chain 042; Cat. #4142) as a substrate in 40 mM sodium acetate, pH 5.5, at 37°C for 20 min. Stop the reaction by adding 1% trifluoroacetic acid. Analyze the reaction mixture by HPLC and calculate the enzyme activity from the amount of the product; PA-lactose (PA-Sugar Chain 026; Cat. #4126).

Properties:

Molecular weight; 60,000 (SDS-PAGE)

Stabilizer; non-ionic detergent (Brij58, NP-40)
bovine serum albumin

Inhibitors; Hg²⁺, Cu²⁺, Fe³⁺, lacto-*N*-biose

Optimum pH; pH 5.5 (Sodium citrate buffer, Sodium phosphate
buffer) Substrate: PA-lacto-*N*-tetraose

pH stability; pH 4.0 - 10.0 (4°C, 16 hours)

Michaelis constant; Km=6.80 μ M
(PA-lacto-*N*-tetraose: PA-Sugar Chain 042, Cat. #4142)
Km=38.9 μ M
(PA-Sugar Chain 003, Cat. #4103)

Purity:

Protease Contamination:

No protease activity is detected after incubation of 10 μ l Lacto-*N*-biosidase with 0.4 mM oxidized insulin B chain for 16 hours at 37°C, in 10 μ l of 200 mM sodium acetate buffer, pH 5.5.

Endoglycosidase and Exoglycosidase Contamination:

Substrate:

1) *p*-NP-glycoside:

25 μ l of Lacto-*N*-biosidase is incubated with 1.1 mM *p*-nitrophenyl glycoside for 16 hours at 37°C, in 225 μ l of 40 mM sodium acetate buffer, pH 4.5. Then 250 μ l of 1 M sodium carbonate is added to the reaction mixture to stop the reaction, measured at 405 nm. The detection limit is 1 μ U/ml.

2) PA-Sugar Chain:

2 μ l of Lacto-*N*-biosidase is incubated with 2 μ M PA-Sugar Chain for 16 hours at 37°C, in 8 μ l of 40 mM sodium acetate buffer, pH 5.5. Then 2 μ l of 1% trifluoroacetic acid is added to the reaction mixture to stop the reaction. The reaction mixture is analyzed by HPLC. The detection limit is 0.1 μ U/ml.

Contaminating activity	<i>p</i> -NP-glycoside	PA-Sugar Chain
α -Fucosidase	< 0.005%	< 0.01%
α -Galactosidase	< 0.02%	ND
β -Galactosidase	ND	ND
α -Glucosidase	NT	ND
β -Glucosidase	< 0.005%	< 0.015%
β - <i>N</i> -Acetylhexosaminidase	< 0.5%	ND
α - <i>N</i> -Acetylgalactosaminidase	NT	ND
α -Mannosidase	ND	ND
β -Mannosidase	< 0.005	NT
Sialidase	NT	ND
Endo- β - <i>N</i> -acetylglucosaminidase	NT	ND

ND; Not detected, NT; Not tested

References:

- 1) Sano M, Hayakawa K, and Kato I. *Proc Natl Acad Sc USA.* (1992) **89**: 8512-8516.
- 2) Sano M, Hayakawa K, and Kato I. *J Biol Chem.* (1993) **268**: 18560-18566.
- 3) Takasaki S and Kobata A. *Biochemistry.* (1986) **25**: 5709-5715.
- 4) Townsend RR, Hardy M R, Wong T C, and Lee Y C. *Biochemistry.* (1986) **25**: 5716-5725.
- 5) Stroud M R, Levery S B, Nudelman E D, Salyan M E K, Towell J A, Roberts C E, Watanabe M, and Hakomori S. *J Biol Chem.* (1991) **266**: 8439-8446.
- 6) Sano M, Hayakawa K, and Kato I. *J Biol Chem.* (1992) **267**: 1522 -1527.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Lacto-*N*-biosidase from *Streptomyces* sp. 142

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Substrate specificity of Lacto-*N*-biosidase

Substrate	Relative activity (%)
Glc β 1-3GlcNAc β 1-3Gal β 1-4Glc-PA (PA-Sugar Chain 042)	100
Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-PA (PA-Sugar Chain 041)	0
Fuc α 1-2Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-PA (PA-Sugar Chain 043)	0
Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-PA Fuc α 1 ⁻⁴ (PA-Sugar Chain 044)	0
Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-PA Fuc α 1 ⁻³ (PA-Sugar Chain 045)	0
Gal β 1-4GlcNAc β 1-2Man α 1 ⁻⁶ Gal β 1-3GlcNAc β 1 ⁻⁴ Man α 1 ⁻³ Man β 1-4GlcNAc β 1-4GlcNAc-PA Gal β 1-4GlcNAc β 1 ⁻² (PA-Sugar Chain 003)	0

Neu5Ac α 2-3 Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-PA 0

PA: Pyridylamino (group)

PA-sugar chain was derivatized according to the method of Kondo, *et al.*

<Application 1 : Digestion of Bovine Asialofetuin by Lacto-*N*-biosidase>

Bovine fetuin is a bovine fetal serum derived glycoprotein (MW48,000), which has *N*-linked and *O*-linked glycosaccharides. It is reported that the structure of *N*-linked glycosaccharide is biantennary or triantennary sugar chain, combined with sialic acid on nonreducing terminus. 20 - 30% of triantennary sugar chain has type I form on nonreducing terminus. The sugar chain structure was analysed with HPLC when asialofetuin, removed sialic acid from fetuin, was treated by Lacto-*N*-biosidase.

[Protocol]

Lacto-*N*-biosidase and 100 pmol asialofetuin were incubated for 2 hours at 37°C in 10 μ l of 40 mM sodium acetate buffer, pH 5.5, then the reaction mixture was boiled to stop the reaction. After evaporation to dryness, 2.5 mU of Glycopeptidase F (Cat. #4450) in 10 μ l of 40 mM sodium phosphate buffer, pH 8.5, was added and incubated for 2 hours at 37°C to release *N*-glycans, then the reaction mixture was boiled to stop the reaction. The reaction mixture was lyophilized, then the released *N*-linked glycosaccharides were analyzed by HPLC. PALPAK Type N (discontinued) and PALPAK Type R (discontinued) were used for analysis.

[Result]

The sugar chain composition of asialofetuin was changed by Lacto-*N*-biosidase treatment. From the analysis result with PALPAK Type R, the rate of PA-Sugar chain 001,002,003 of non-treated asialofetuin was 1 : 6 : 3. But the rate of treated asialofetuin was changed to 3 : 6 : 1. It showed that almost triantennary sugar chains with the type I chain were changed into biantennary sugar chains.

<Application 2 : Digestion of Le^a/Le^a glycolipid by Lacto-*N*-biosidase>

Le^a/Le^a glycolipid is a kind of lact series glycolipid which structure contains type I sugar chain repeat. This glycolipid was digested with endoglycoceramidase to release sugar chains, and it was pyridylaminated by *GlycoTAG*. This pyridylaminated sugar chain was digested by α -1,3/4-L-Fucosidase (Cat. #4453) and Lacto-*N*-biosidase. The product was analyzed by HPLC using PALPAK Type N.


[Protocol]

All reaction was performed at 37°C for 30 min in 30 mM sodium acetate, pH 5.5. Amount of enzyme for each reaction was 0.05 μ U/pmol sugar chain of α -1,3/4-L-Fucosidase or 0.5 μ U/pmol sugar chain of Lacto-*N*-biosidase.

- 1) PA-Le^a/Le^a was digested with α -1,3/4-L-Fucosidase, then boiled to stop the reaction. 10 pmol was analyzed by HPLC.
- 2) Lacto-*N*-biosidase was added to the resting solution and performed, then reaction mixture was boiled to stop the reaction. 10 pmol was analyzed by HPLC.
- 3) α -1,3/4-L-Fucosidase was added to the resting solution and performed, then the reaction mixture was boiled to stop the reaction. 10 pmol was analyzed by HPLC.
- 4) Lacto-*N*-biosidase was added to the resting solution and performed, then the reaction mixture was boiled to stop the reaction. 10 pmol was analyzed by HPLC.

[Result]

The result of HPLC analysis showed that PA-Le^a/Le^a was sequential digested with α -1,3/4-L-Fucosidase and Lacto-*N*-biosidase.

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