qPCR Subtyping Primer Set for NA of Avian Influenza

Code No. 6641

Size: 100 μ l each

Conc.: 5X

Description :

This product is a primer set, which is designed based on nucleotide sequence of Avian influenza virus NA (Neuraminidase) gene. Using this set with PrimeScript™ RT reagent Kit (Perfect Real Time) (Cat. #RR037A/B) and TB Green® *Premix Ex Taq*™ (Tli RNaseH Plus) (Cat. #RR420A/B)*¹ allows analytical determination for NA subtypes of Avian influenza virus by RT-qPCR.*²

- * 1 We have begun the process of changing the names for Takara Bio's intercalator-based real-time PCR (qPCR) products to the "TB Green series". These products can be used the same way as before, as only the names are changing. Catalog number and product performance are unaffected by this transition.
- * 2 There might be genes that could be difficult for this analytical determination, although this product is optimized for many NA subtypes.

Takara Bio received a license from the National Agriculture and Food Research Organization, and commercialized this product.

Storage: -20°C

Subtypes for detection: N1 ~ N9

Target gene :

- 1 - 11. NA (Neuraminidase) gene
- 12. NP (Nucleoprotein) gene

Contents (for 25 reactions in each primer) :

Following 12 primer solutions (Forward & Reverse primer mix) are included. (100 μ l each)

Form: Solution in TE Buffer (10 mM Tris-HCl (pH 8.0), 0.1 mM EDTA)

1. N1 Subtype Eu/Am Lineage Primer
2. N2 Subtype Eu/Am Lineage Primer
3. N3 Subtype Eu/Am Lineage Primer
4. N4 Subtype Eu/Am Lineage Primer
5. N5 Subtype Eu/Am Lineage Primer
6. N6 Subtype Eu/Am Lineage Primer
7. N7 Subtype Eu/Am Lineage Primer
8. N8 Subtype Eu Lineage Primer
9. N8 Subtype Am Lineage Primer
10. N9 Subtype Eu Lineage Primer
11. N9 Subtype Am Lineage Primer
12. NP Eu/Am Lineage Primer

Eu: Eurasian lineage
Am: North American lineage
Eu/Am: Eurasian lineage and North American lineage

Precautions before use:

This product is designed for detection of the target gene DNA, and detects not only living virus but also inactivated virus. The target gene may not be detected when there is mutation, deletion, or insertion in the genomic sequences correspond with the primers. (Takara Bio is not responsible for any actions taken as a result of analytical determinations made with this product.)

Application example :

< Reverse transcription reaction >

- Reagent : PrimeScript RT reagent Kit (Perfect Real Time)
- Reaction mixture components :

5X PrimeScript Buffer (for Real Time)	2 μ l
PrimeScript RT Enzyme Mix I	0.5 μ l
Random 6 mers (100 μ M)	0.5 μ l
total RNA	
RNase Free dH ₂ O	
Total	10 μ l* ³

- * 3 500 ng of total RNA is the maximum amount for the reverse transcription in 10 μ l reaction mixture. It is possible to scale up the reverse transcription reaction as required.

- Reaction condition for reverse transcription :

37°C	15 min
85°C	5 sec
4°C	

< qPCR reaction >

- Reagent : TB Green *Premix Ex Taq* (Tli RNaseH Plus)
- Reaction mixture components :

The reaction mixture has to be prepared for each primer solution.

TB Green <i>Premix Ex Taq</i> (Tli RNaseH Plus)	10 μ l
Primer	4 μ l
RT reaction solution (cDNA)	1 μ l
Sterile purified water	5 μ l
Total	20 μ l

- qPCR condition :

95°C	30 sec	} 30 cycles
95°C	30 sec	
50°C	20 sec	
72°C	32 sec	

Dissociation curve analysis

Note : Addition of ROX Reference Dye is required in reaction mixture when using real-time PCR instruments of Applied Biosystems. Please see the instruction for TB Green *Premix Ex Taq* (Tli RNaseH Plus) for more detail.

Examples for reaction result :

Reaction results (example) using each primer are shown on Takara Bio web catalog.

Reference :

Tsakamoto K, Javier PC, Shishido M, Noguchi D, *et al.* *J Clin Microbiol.* (2012) **50**: 37-45.

TB Green is a registered trademark of Takara Bio Inc.
PrimeScript and *Premix Ex Taq* are trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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