Anhydrotrypsin Agarose

Code No. 7302

Size: 1 ml wet gel

Form :

Gel suspension in 50 mM sodium acetate buffer, pH 5.0, containing 20 mM CaCl₂ and 0.02% NaN₃

Description : Anhydrotrypsin Agarose is an affinity chromatography adsorbent that selectively adsorbs peptides with Arg, Lys, or AECys (S-aminoethyl cysteine) at the C-termini under mildly acidic conditions. It does not bind free amino acids, or peptides which have Arg, Lys, or AECys at sites other than the C-termini.

By use of this specific adsorbent, it is possible to isolate C-terminal peptide fragments from a protease digest of proteins.

Storage: 4°C

Do not keep the gel in elution buffer (such as HCl or formic acid solutions, pH 2.5) for extended periods.

Binding Capacity:

Soybean trypsin inhibitor binding capacity:

approx. 80 nmol/ml wet gel

Bz-Gly- Arg binding capacity:

approx. 60 nmol/ml wet gel

Note :

This product contains sodium azide as a preservative. Sodium azide reacts with many heavy metals to form explosive compounds. It is critical that a large amount of water be used when disposing.

References:

- 1) Ishii S, Yokosawa H, Kumazaki T, and Nakamura I. *Methods in Enzymology*. (1983) **91**: 378.
- 2) Kumazaki T, Nakako T, Arisaka F, and Ishii S. *PROTEINS: Structure, Function, and Genetics*. (1986) **1**: 100.
- 3) Yokosawa H, and Ishii S. *Anal Biochem*. (1979) **98**: 198.
- 4) Arisaka F, Nakako T, Kumazaki T, and Ishii S. *J Protein Chem*. (1987) **6**: 245.
- 5) Kumazaki T, Terasawa K, and Ishii S. *J Biochem*. (1987) **102**: 1539.

General Method:

1. Preparation of Sample Prior to Chromatography

- A protein sample is digested by trypsin (or by chymotrypsin). The reaction is terminated by adding DFP (diisopropyl fluorophosphate, final conc. 1 mM), and the pH of the digest is adjusted to 5.0 . No need to add 0.02 M CaCl₂
- The sample is kept to a smaller volume than the gel volume.
- The amount of protease digest to be applied to the agarose column should be determined by referring to the Bz-Gly-Arg binding capacity.
- Samples containing 2 - 4 M urea could be directly applied to the agarose column after adjusting to pH 5.0.

2. Method for chromatography

The following procedures should be done at 4°C.

1. The gel is put into a column, washed with 0.05 M sodium acetate buffer, pH 5.0, containing 0.02 M CaCl₂ (starting buffer). Washing volume should be over 20 - fold gel bed volume.
2. The sample is applied onto this column.
3. The column is washed with starting buffer (20 - fold of gel bed volume) at flow rate of about 17 ml/hr/cm². Peptides with Arg, Lys, or AECys C-terminal residues are adsorbed. Peptides with Arg-Pro or Lys-Pro at C-termini are not adsorbed though trypsin does not cleave such peptide bonds.
4. Peptides bound to the gel can be eluted by 5 mM HCl (approx. pH 2.5) or 0.1 M HCOOH (approx. pH 2.5). Sufficient recovery is achieved by elution with HCl or HCOOH at volumes of 10 - fold that of the gel.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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