# Calpastatin

Code No. 7316

Size:

3 mg

## Description :

Calpastatin is an endogenous protease inhibitor that acts specifically on calpain (a calcium-dependent cysteine protease). It contains four repetitive sequences composing of 120 to 140 amino acid residues, and an N-terminal non-homologous sequence.<sup>1,2)</sup>

This product is highly purified recombinant human calpastatin domain 1, which consists of 135 amino acids (MW 14,000).<sup>1,3)</sup>

**Storage:** -20°C

**Form:** Lyophilized white powder  
(Reconstitute with sterile purified water.)

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Typical Activity :

50 nM of calpastatin domain 1 completely inhibits the activity of 7.5 µg/ml calpain I. 15 nM of calpastatin domain 1 inhibits 50% of the activity of 7.5 µg/ml calpain I.

## References :

- 1) Asada K, Ishino Y, Shimada M, Shimojo T, Endo M, Kimizuka F, Kato I, Maki M, Hatanaka M, and Murachi T. *J Enzyme Inhibition*. (1989) **3**: 49-56.
- 2) Kanaki R and Murachi T. *Protein, nucleic acid and enzyme*. (1987) **32**: 116-129 (Japanese Journal).
- 3) Uemori T, Shimojo T, Asada K, Asano T, Kimizuka F, Kato I, Maki M, Hatanaka M, Murachi T, Hanzaw H, and Arata Y. *Biochem Biophys Res Comm*. (1990) **166**: 1485-1493.

## Assay of activity :

### Principle

The activity of calpastatin is analyzed as an inhibition activity against calpain.

### Protocol

(Reagents)

- A. 1 mg/ml Calpain
- B. 2% Casein
- C. 6 mg/ml Cysteine
- D. 0.03 - 0.05 mg/ml Calpastatin in H<sub>2</sub>O
- E. 50 mM CaCl<sub>2</sub>
- F. 5% TCA

(Procedure)

1. Mix 1.5 µl of Reagent A, 40 µl of Reagent B, 20 µl of Reagent C, and 10 µl of Reagent D.
2. Incubate for 5 min at 37°C.
3. Add 20 µl of Reagent E and incubate at 37°C for 20 min.
4. Terminate the reaction by adding 200 µl of Reagent F.
5. Centrifuge, and read absorbance of the supernatant at 280 nm.
6. For controls, add H<sub>2</sub>O instead of Reagent D (Blank 1) and also for Reagent D and E (Blank 2).

(Calculations)

Inhibition activity (µg calpain inhibited/µg Calpastatin)

$$= \frac{\Delta A_{280 \text{ Blank}} - \Delta A_{280 \text{ Test}}^*}{\Delta A_{280 \text{ Blank}}} \times \frac{\text{Calpain } (\mu\text{g})}{\text{Calpastatin } (\mu\text{g})}$$

$$\Delta A_{280 \text{ Blank}} = A_{280 \text{ Blank 1}} - A_{280 \text{ Blank 2}}$$

$$\Delta A_{280 \text{ Test}} = A_{280 \text{ Test}} - A_{280 \text{ Blank 2}}$$

\* Adjust the concentration of calpastatin to give  $\Delta A_{280 \text{ Test}} = 0.1 - 0.3$ .

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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