# Endoproteinase Asp-N

**Code No. 7329**      **Size:**      **2  $\mu$ g**

**Supplied Reagent (5X):**  
**250 mM Sodium phosphate, pH 8.0**      **1 ml**

## Description:

Endoproteinase Asp-N is a metalloprotease which specifically cleaves the peptide bonds at the amino side of aspartic acid or cysteine residues of proteins and peptides.

**Source:** *Pseudomonas fragi* mutant

**Form:** Lyophilized  
(containing the equivalent to 50  $\mu$ l of 10 mM Tris-HCl, pH 7.5)  
Reconstitute lyophilized Endoproteinase Asp-N in 50  $\mu$ l of distilled water.

**Storage:** below 4°C (under dry condition)

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Definition of Activity:

One azocoll unit is defined as the activity that will set free as much azo dyestuff which yields an absorbance of 0.001 at 520 nm per minute at 37°C, pH 7.5.

**Specific Activity:**  $\geq 20$  U/ $\mu$ g

## Properties:

Optimum pH: 7.0 - 8.0

Optimum temp: 37°C

Inhibitors: 2-phenanthroline, EDTA

Tolerance to denaturants:

- Stable against  $\leq 1$  M urea
- $\leq 1$  M guanidine-HCl
- $\leq 0.01\%$  SDS
- $\leq 10\%$  (v/v) acetonitrile

## References:

- 1) Noreau J and Drapeau G R. *J Bacteriol.* (1979) **140**: 911-916.
- 2) Drapeau G R. *J Biol Chem.* (1980) **255**: 839-840.
- 3) Geuss U, et al. *J Protein Chem.* (1990) **9**: 299-300.

## Application:

Cleavage of reduced soluble lysozyme by Endoproteinase Asp-N

## Reaction Condition

Substrate: 1 nmol  
Enzyme: 0.25  $\mu$ g (approximately 10 pmol)  
Reaction buffer: 50 mM sodium phosphate, pH 8.0  
Reaction volume: 100  $\mu$ l  
Reaction temperature: 37°C  
Reaction time: 4 hours

## Results:

Results are shown in Figure 1 and 2. Non-specific cleavage at sites other than -X-Asp- was not observed.


Figure 1. HPLC of digests of reduced soluble lysozyme by Endoproteinase Asp-N

Column: Nucleosil 5C<sub>18</sub>-AR (4.6  $\phi$  x 150 mm)  
Solvent A: 0.05% TFA/H<sub>2</sub>O  
Solvent B: 0.05% TFA/CH<sub>3</sub>CN  
Elution: 0  $\rightarrow$  15% B/0  $\rightarrow$  15 min  
15  $\rightarrow$  27.5% B/15  $\rightarrow$  65 min  
27.5  $\rightarrow$  42.5% B/65  $\rightarrow$  80 min  
Flow rate: 0.5 ml/min  
Detection: 215 nm


Figure 2. Amino acid sequence of cleaved peptide fragments

Peptide fragments produced from cleavage by Endoproteinase Asp-N are shown. The numbers under the line correspond to the peaks in HPLC (Figure 1).

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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