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The screenshot shows the Exiqon website header with the logo 'EXIQON Seek Find Verify' and navigation links: Home, Search, Contact, Print, Basket, and Sign In. The Sign In link is circled in red and labeled with a circled '1'. Below the header is a navigation menu with links: Products, Services, Applications, Resource Center, Ordering, About Exiqon, and Investor. The main content area is titled 'Sign in' and contains a form with the following elements: a 'Sign in' logo, 'Username or email:' and 'Password:' input fields, a 'Remember me next time.' checkbox, a 'Sign in' button, and 'Register' and 'Forgot password?' links. The 'Register' link is circled in red and labeled with a circled '2'.

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[Home](#) > [microRNA PCR](#) > Pick-&Mix microRNA PCR Panels

Pick-&Mix microRNA PCR Panels

Design your own 96- and 384-well qPCR plates. Drag and drop miRCURY LNA™ Universal RT microRNA PCR assays or any custom microRNA, mRNA or lncRNA assay onto the plates using our easy-to-use online design tool. Pick-&Mix the assays most relevant for your project.

[Configure your plates](#)

IMPORTANT NOTICE

As of September 16 2014, Exiqon converted to manufacture Pick-&Mix panels based on the V4 based assay lists. This may have potential implications on your ongoing microRNA studies using Exiqon's V3 products. If you have already saved Pick-&Mix panel configurations in your Exiqon Profile or by other means configured panel designs with Exiqon's Pick-&Mix assay lists previous to this date, you are requested to check if some obsolete V3 assays may have been removed from your layout and if that is the case, redesign your panel. Please contact us if you have any questions.

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[Exit configurator](#)

[View Saved Panels](#)

Welcome to the Pick&Mix plate configuration tool

Step 1

During the next 7 steps you will be able to configure your Pick&Mix plates by selecting the appropriate plate type for your qPCR instrument, target organism, plate layout, as well as selecting and positioning the LNA™ microRNA PCR primer sets of choice.

- ① Please select or confirm your desired plate format below (please change your selection if this is not compatible with your qPCR instrument).

96
 384

The Pick&Mix plates are priced according to plate format (96 well or 384 well plates) regardless of which layout or which microRNA primer sets you choose. The minimum batch size is 8 plates with the option of adding additional plates to the order in increments of 4 plates. For prices, please click [here](#).

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The screenshot shows the 'New Plate' configuration interface. At the top, there is a blue header bar with the text 'New Plate'. Below the header, there is a progress bar with six steps: Step 1/6, Step 2/6 (selected), Step 3/6, Step 4/6, Step 5/6, and Step 6/6. To the right of the progress bar are two buttons: 'Exit configurator' and 'View Saved Panels'. Below the progress bar, the section is titled 'Step 2' and 'Select qPCR instrument'. The text below reads: 'The Pick&Mix product is available in various plate types to provide compatibility with a wide range of qPCR instruments. Select the instrument for which your Pick&Mix plates need to be compatible with. If your qPCR instrument does not appear in the list, please [contact us](#).' There are four radio button options: 'ABI 7900HT Real-Time PCR system' (selected), 'ABI ViiA 7 (384-well block) Real-Time PCR system', 'Bio-Rad CFX 384 Real-Time PCR system', and 'Roche LightCycler 480 Real-Time PCR system'. At the bottom, there are two buttons: 'Previous' and 'Next' (selected).

New Plate

Step 1/6 Step 2/6 Step 3/6 Step 4/6 Step 5/6 Step 6/6 Exit configurator View Saved Panels

Step 2

Select qPCR instrument

The Pick&Mix product is available in various plate types to provide compatibility with a wide range of qPCR instruments. Select the instrument for which your Pick&Mix plates need to be compatible with. If your qPCR instrument does not appear in the list, please [contact us](#).

①

- ABI 7900HT Real-Time PCR system
- ABI ViiA 7 (384-well block) Real-Time PCR system
- Bio-Rad CFX 384 Real-Time PCR system
- Roche LightCycler 480 Real-Time PCR system

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New Plate

▶ Step 1/6 ▶ Step 2/6 **▶ Step 3/6** ▶ Step 4/6 ▶ Step 5/6 ▶ Step 6/6 [Exit configurator](#) [View Saved Panels](#)

Step 3

Select target organism and assay list
Select the type of organism and then which list of LNA™ microRNA primer sets you would like to use in designing your panel. You may select the full list of assays available for your organism of interest. Alternatively, choose list of selected assays based on one of Exiqon's pre-designed panels.

① **Human** Mouse Rat More Organisms

Select assay list:

② All assays

Assay list of pre-designed panels:

- Human Panel I (miRNome, 372 human microRNAs, 3 ncRNA controls) [details](#) ▶
- Cancer Focus Panel (84 human microRNAs, 3 ncRNA controls) [details](#) ▶
- Serum Plasma Focus Panel (179 human microRNAs) [details](#) ▶
- Stem Cell (Pluripotent) Focus Panel (85 microRNAs, 3 nc RNA controls) [details](#) ▶
- Toxicology Focus Panel [details](#) ▶
- Viral assays [details](#) ▶
- Urine Exosome Focus Panel [details](#) ▶
- CSF Exosome Focus Panel [details](#) ▶

For other organisms please [contact us](#). Alternatively, you may view the [list of pre-designed microRNA primer sets](#) for determining which set of assays will best cover your needs.

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New Plate

Step 1/6 Step 2/6 Step 3/6 **Step 4/6** Step 5/6 Step 6/6

Exit configurator View Saved Panels

Step 4

Name panel and select plate layout
Select a plate layout and assign a unique name to your plate. [More details](#)

In order for us to provide you with the final plate layout electronically, we take your email address from logged user profile. You can provide alternative email address as well. [More details](#)

Please provide a name for the plate: (①)

Alternative email address (optional):

Select fully flexible or pre-designed plate layout
Select one of the different plate layouts depicted below. You may choose to design in a fully flexible plate layout placing your assays in a pattern and with replicates of your wish. Alternatively, you may select one of the pre-defined plate layouts shown to the right for easy assay replicate positioning. Click on the plate layouts for enlarged view. [More details](#)

②

- 380 microRNAs x 1 sample
 - 1 x 4 UniSp3 IPC
 - 1 x 16 UniSp6 CP (optional)
- 22 microRNAs x 16 samples
 - 1 x 16 UniSp3 IPC
 - 1 x 16 UniSp6 CP (optional)
- 46 microRNAs x 8 samples
 - 1 x 8 UniSp3 IPC
 - 1 x 8 UniSp6 CP (optional)
- 94 microRNAs x 4 samples
 - 1 x 4 UniSp3 IPC
 - 1 x 4 UniSp6 CP (optional)

Ideal layout for 96-well multi-channel pipette

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① primer set の検索

Select from: Search:

Assign all pre-defined assays to panel:

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▶ Step 1/6 ▶ Step 2/6 ▶ Step 3/6 ▶ Step 4/6 ▶ Step 5/6 Step 6/6

Exit configurator View Saved Panels Configure new Panel

Step 6

Download plate layout. Confirm and finalize panel configuration
Your panel has been saved. Please download your draft plate layout and review the layout before proceeding.

Your plate layout is downloaded as an Excel file. For easy identification, it is recommended to keep the file name when downloading. The file is designed for use in your experimental setup and subsequent data analysis.

The Plate layout file also serves as your GenEx template file. It is fully compatible with Exiqon's GenEx qPCR analysis software for direct import into the software in the present format for easy normalization and data analysis. [Read More](#)

Proceed by selecting one of the following three options:

.....

Save panel and exit panel configurator. Click save panel if you are done configuring the Panel for now and exit the panel configurator without further handling. You may at any time retrieve you panel from your personal My Profile account if you wish to edit or order the panel.

.....

Make panel corrections. Click edit panel to make corrections to the panel if you need to make changes to the instrument, plate or assay selection for the current panel before saving or ordering the panel.


.....

Ready to order panel. Click finalize panel and proceed to check-out if you are satisfied with the panel and would like to order or prepare for ordering the panel. The panel design will be locked from further editing. The next step will allow you to:

- > Select number of plates and see final pricing
- > Ordering directly online or by downloading a web shop PO, required if contacting Exiqon for a quote or when placing an order by PO, email, or fax

An email with the plate layout will be sent to the email address provided in your personal My Profile account. Please save this plate layout or download it at any time from your My Profile account. You will not receive a copy of the plate layout with the delivery of your Pick-&-Mix microRNA PCR Panel.

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Basket

Since you are located in a distributor country, Exiqon cannot offer you price information. To receive a quote, simply press the "Get quote" button below. You will be taken to a new page where you will be asked to submit your contact details. Your local distributor will then contact you with a quote.

TEST layout

Product number: 203819

Description: miRCURY LNA™ Universal RT microRNA PCR, Ready-to-use PCR panels with custom selection of LNA™ microRNA PCR primer sets, 384-well PCR Plate First 8 plates. For use with ABI 7900HT, Vii7 and Bio-Rad CFX384., Instrument: ABI 7900HT Real-Time PCR system, Layout: 380 microRNAs x 1 sample

Plate No: 675851375

Remove

Quantity 8 plates

The following products are required with the items in your shopping basket:

Product no.: 203301

Product name: Universal cDNA Synthesis Kit II, 8-64 rxns

Description: miRCURY LNA™ microRNA PCR, Polyadenylation and cDNA synthesis kit II (8-64 rxns)

Please inquire

Product no.: 203402

Product name: ExiLent SYBR® Green master mix, 2.5ml

Description: miRCURY LNA™ microRNA PCR, ExiLent SYBR® Green master mix, 2.5 ml (500 rxns of 10ul)

Please inquire

Product no.: 203420

Product name: ExiLent SYBR® Green master mix, 20ml

Description: miRCURY LNA™ microRNA PCR, ExiLent SYBR® Green master mix, 20 ml (4000 rxns of 10ul)

Please inquire

Press the 'Check out' button to continue the check out procedure online.

If you would rather like to email/fax your order, please follow these instructions:

→ Press the 'Save quote - Order later' button.

→ Email/Fax us this quote along with the following information:

- Your name
- Billing address
- Shipping address
- A copy of the purchase order
- TAX/VAT ID number

→ Email to: order@exiqon.com or use the following fax numbers:

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Contact us

In our effort to give you a rapid and fulfilling answer, we kindly ask you to fill in a few blanks below

Your inquiry is related to

- | | |
|--|---|
| <input type="radio"/> Technical Support | <input type="radio"/> Career in Exiqon |
| <input checked="" type="radio"/> Product Information | <input type="radio"/> Web & Marketing |
| <input type="radio"/> Ordering | <input type="radio"/> Sign up to Newsletter |
| <input type="radio"/> Profiling Services | <input type="radio"/> Other |

- Sign up to our Newsletter and get research news, free downloads and offers 8-10 times a year. You can unsubscribe at anytime.

Please chose area of interest (click all that apply)

- | | |
|---|--|
| <input checked="" type="checkbox"/> microRNA Research | <input type="checkbox"/> mRNA Research |
| <input type="checkbox"/> Molecular Diagnostics | <input type="checkbox"/> Other |

Please provide us with information about your inquiry. We will return to you as soon as the inquiry has been processed

Please send me a product quote. Reference number:

 [Click here to attach a file](#)

Reference number

xxxxxx

Please state contact details

First name	<input type="text"/>	*
Last name	<input type="text"/>	*
E-mail	<input type="text"/>	*
Phone	<input type="text"/>	
Institution	<input type="text"/>	*
Address	<input type="text"/>	*
City	<input type="text"/>	*
Zip	<input type="text"/>	*
Country		*
In which industry are you employed		*
Are you currently engaged in microRNA research		*
How did you first learn about Exiqon		*

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