Westase

(Digestive enzyme of cell walls of yeast)

Code No. 9005

Size : 1 g

Shipping at – 20°C
Store at 4°C

Lot No.

Expiration Date :

Quality :

β -1,6 glucanase activity : U/g powder
Lytic activity : U/g powder

Origin : *Streptomyces rochei* DB-34

Form : Lyophilized powder (containing celite as the excipient)

Storage : 4°C , dry condition.

Description :

This product was prepared from liquid culture supernatant of *Streptomyces rochei* DB-34. This product has complex lytic activities of yeast cell mainly consisting of β -1,6 glucanase and β -1,3 glucanase activity.

Unit definition :

β -1,6 glucanase activity : One unit is defined as the amount of enzyme required to release 1 μ mol reducing sugar from 10 mg/ml Pustulan solution in 1 min. at 37°C , pH6.0.

Lytic activity : One unit is defined as the amount of enzyme required to cause a 1% decrease in absorbance at 660 nm in 1 min. at 30°C, pH6.0 when using cell wall fraction of *Cryptococcus albidus* IFO 0612 as a substrate.

Specifications :

β -1,6 glucanase activity : \geq 400 units/g
Lytic activity : \geq 35,000 units/g
DNase activity : ND (McIlvain Buffer, pH6.0)

Preparation of the enzyme solution :

- (1) Dissolve in McIlvain Buffer*.
*: McIlvain Buffer is prepared by mixing 0.1M Citric acid solution and 0.2 M Disodium hydrogenphosphate at the ratio of around 36.8 : 63.2 (v/v), and adjust to pH6.0.
- (2) Filter the solution with cellulose acetate filter.

Note:

- (1) Na-tartrate (0.3 - 0.5 M) must be used as an osmotic stabilizer to form protoplast by Westase. The use of sorbitol as a stabilizer remarkably reduces the efficiency of forming protoplast.
- (2) For preparation of protoplasts by Westase, the yeast cells in exponential phase are suitable. In most yeast strains, the cells of stationary phase are not appropriate since they can result in lower efficiency.
- (3) To regenerate the cell wall of protoplasts, use 1.0 - 1.5 M sorbitol as an osmotic stabilizer.

Example of protoplast preparation :

- (1) *Saccharomyces cerevisiae*, *Trigonopsis variabilis*
Medium YPG medium
Culture 25°C , 1day, Shake culture (reciprocate at 120 rpm)
Conditions 0.5% Westase solution
0.5 M Na-tartrate
McIlvain Buffer (pH6.0)
30°C Reciprocal shaking, 1 - 2 hr.
- (2) *Lipomyces starkeyi*, *Candida utilis*
Medium Same as in (1)
Culture Same as in (1)
Conditions 0.5% Westase solution
0.4 M Na-tartrate
McIlvain Buffer (pH6.0)
30°C Reciprocal shaking, 1 - 2 hr.
- (3) *Shizosaccharomyces pombe*, *Hansenula mrakii*, *Kluyveromyces lactis*, *Pichia anomala*, *Filobasidium floriforme*, *Candida colliculosa*
Medium Same as in (1)
Culture Same as in (1)
Conditions 0.5% Westase solution
0.4 M Na-tartrate
McIlvain Buffer (pH6.0)
30°C Reciprocal shaking, 3 - 4 hr.
- (4) *Kloeckera apiculata*, *Cryptococcus albidus*
Medium Same as in (1)
Culture Same as in (1)
Conditions 0.5% Westase solution
0.3 M Na-tartrate
McIlvain Buffer (pH6.0)
30°C Reciprocal shaking, 1 - 2 hr.
- (5) *Ustilago maydis*, *Graphiophora phoenicis*, *Brettanomyces bruxellensis*, *Phaffia rhodozyma*
Medium Same as in (1)
Culture 25°C , 2days, Shake culture (reciprocate at 120 rpm)
Conditions Same as in (1)
- (6) *Tremella mesenterica*
Medium Same as in (1)
Culture Same as in (5)
Conditions Same as in (2)

Manufactured by Ozeki corporation.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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- (1) *Saccharomyces cerevisiae*, *Trigonopsis variabilis*
培地 YPG medium
培養条件 25°C、1day、Shake culture (reciprocate at 120 rpm)
調製条件 0.5% Westase solution
0.5 M Na-tartrate
Mcllvain Buffer (pH6.0)
30°C Reciprocal shaking、1 - 2 hr.
- (2) *Lipomyces starkeyi*, *Candida utilis*
培地 Same as in (1)
培養条件 Same as in (1)
調製条件 0.5% Westase solution
0.4 M Na-tartrate
Mcllvain Buffer (pH6.0)
30°C Reciprocal shaking、1 - 2 hr.
- (3) *Shizosaccharomyces pombe*, *Hansenula mrakii*, *Kluyveromyces lactis*, *Pichia anomala*, *Filobasidium floriforme*, *Candida colliculosa*
培地 Same as in (1)
培養条件 Same as in (1)
調製条件 0.5% Westase solution
0.4 M Na-tartrate
Mcllvain Buffer (pH6.0)
30°C Reciprocal shaking、3 - 4 hr.
- (4) *Kloeckera apiculata*, *Cryptococcus albidus*
培地 Same as in (1)
培養条件 Same as in (1)
調製条件 0.5% Westase solution
0.3 M Na-tartrate
Mcllvain Buffer (pH6.0)
30°C Reciprocal shaking、1 - 2 hr.
- (5) *Ustilago maydis*, *Graphioli phoenicis*, *Brettanomyces bruxellensis*, *Phaffia rhodozyma*
培地 Same as in (1)
培養条件 25°C、2days、Shake culture (reciprocate at 120 rpm)
調製条件 Same as in (1)
- (6) *Tremella mesenterica*
培地 Same as in (1)
培養条件 Same as in (5)
調製条件 Same as in (2)

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