Proteinase K

Code No. 9034
Size : 5 ml

Shipping at 4°C
Store at 4°C

Lot No.
Expiration Date :

Description :

Proteinase K is a serine protease with high activity. Its enzyme activity is stable under the existence of calcium salt, and increases when supplied with protein denaturant such as SDS or urea.

Proteinase K can be available for a wide range of substrate. Especially it degrades estel and peptide band preferentially next to C-termini of hydrophobic, sulfuric, or aromatic amino acid.

Source :

Tritirachium album

Unit definition :

One Anson U of enzyme activity is defined as the amount required to liberate 1.0 μ mol tyrosine from hemoglobin per minute at 37°C, pH7.5.

Concentration :

> = 600 mAnson U/ml

Purity :

Nuclease activity is not detected in any of the following cases:

1. After incubation of pBR322 with this enzyme for 6 hours at 37°C.
2. After incubation of RNA with this enzyme for 2 hours at 37°C.

Applications :

1. Purification of DNA, RNA, phage
2. Preparation of chromosomal DNA for *in situ* hybridization, fingerprinting, colony hybridization, plaque hybridization and pulsed-field electrophoresis.

Application example :

0.01 M Tris-HCl, pH7.8
0.01 M EDTA
0.5% SDS
1/400 - 1/40 (v/v) Proteinase K

↓

37 - 56°C, 1hr - overnight incubation

↓

Use the mixture for the next reaction, following phenol/chloroform extraction and ethanol precipitation.

References :

- 1) Ebeling, W., *et al.* (1974) *Eur. J. Biochem.*, **47**, 91-97.
- 2) Sambrook, J., *et al.* (1989) *Molecular Cloning: A Laboratory Manual*, 2nd edition, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- 3) Gunkel, F. A. and Gassen, H. G. (1989) *Eur. J. Biochem.*, **179**, 185-194.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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