EASY Dilution (for Real Time PCR)

Code No. 9160

Size:

1 ml x 8

* 2 years from date of receipt under proper storage conditions.

Description:

This product is used to dilute template DNA or RNA when preparing the serial dilution of standard for establishing the standard curve in real time PCR. When template DNA and RNA are diluted with sterile purified water or TE buffer, they can stick to the inside tube which may prevent correct dilution. This product helps to avoid occurring this absorption and the correct dilution at lower concentrations is possible. In addition, as this product does not include nucleic acid like tRNA or rRNA, non-specific amplification derived from those RNA will not happen.

Storage: -20°C

Note:

EASY Dilution is designed to use with real-time PCR reagents produced by Takara Bio. The combination of this product and other company's reagents hasn't been tested.

Protocol:

The concentration of this product is 1X. It can be directly used for dilution of template DNA or RNA. Serial dilution of the standard solution is added into reverse transcription or PCR reaction at an arbitrary amount. But it is recommended to add a certain volume in a series of experiments.

(Example)

Preparation for serial dilution of cDNA* for drawing the standard curve:

Prepare the serial dilution series of cDNA at the concentration of 5 ng/μl, 500 pg/μl, 50 pg/μl, 5 pg/μl, and 0.5 pg/μl by diluting a portion of the cDNA solution (50 ng/μl) by 10-fold each.

* The cDNA was synthesized from total RNA through reverse transcription.

By adding 2 μl of serial dilution to each reaction, the standard curve toward the cDNA in the range of 1 pg to 100 ng is drawn. (For negative control reaction, add this product by 2 μl each.)

Note:

It is recommended to store this product in dispensing into aliquots convenient for each customer in order to avoid contamination among samples.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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