High-Salt Solution for Precipitation (Plant)

Code No. 9193

Size:

50 ml

* 2 years from date of receipt under proper storage conditions when unused. Once opened, it should be used earlier avoiding contamination.

Description :

High-Salt Solution for Precipitation (Plant) is a reagent for removal of polysaccharides from RNA solution which extracted from plants containing a lot of polysaccharides and proteoglycans, *et al.* High purity of RNA sample is easily obtained only by addition of this product at isopropanol precipitation in purification steps of RNA extraction procedure.

Form : 1.2 M NaCl
0.8 M Sodium Citrate

Storage: Room Temperature

Usage : Removal of polysaccharides from RNA solution.

Purity : DNase and RNase were not detected in this product.

Protocol :

- 1) Add 0.5 volume of this solution to RNA sample*¹, and mix.
- 2) Add same volume of isopropanol as the solution added at 1), and mix.
- 3) Stand for 10 minutes at room temperature.
- 4) Centrifuge at 12,000 rpm for 10 minutes at 4°C.
- 5) Discard the supernatant.
- 6) Add 1 ml of cold 75% ethanol to the precipitate.
- 7) Re-centrifuge at 12,000 rpm for 5 minutes at 4°C.
- 8) Discard the supernatant, and dry the precipitate briefly.*²
- 9) Dissolve the precipitate in appropriate amount of RNase-free water or TE buffer.

* 1 This product is also available for the extraction solution before isopropanol precipitation during RNA extraction procedure by using RNAiso Plus (Cat. #9108/9109).

* 2 Do not dry completely the precipitate. It may cause difficulty with dissolving RNA.


Notes :

1. It is not recommended to use for recovery of small molecular RNA, because of low recovery of that.
2. When RNA sample is small amount, it is recommended to use Dr. GenTLE™ Precipitation Carrier (Cat. #9094) at isopropanol precipitation for increasing the RNA recovery.

Application example :

Three samples of RNA solutions prepared from tobacco cultured cell (BY2), of which A_{260}/A_{230} ratio are less than 1.5, are treated with this product for removal of polysaccharides according to the Protocol. After treatment, the A_{260}/A_{230} ratio of the RNA samples became more than 2.0, showing that the polysaccharides in the RNA samples are removed.

Removal of polysaccharides

Related product :

Dr. GenTLE™ Precipitation Carrier (Cat. #9094)
RNAiso Plus (Cat. #9108/9109)

Dr. GenTLE is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

v201902Da

High-Salt Solution for Precipitation (Plant)

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