TaKaRa Taq™ HS Low DNA

Code No. R090A Size: 200 µl x 5
(for 100 PCR reactions)

Supplied Reagent:
dH₂O for Low DNA 200 µl x 5

Description :

TaKaRa Taq HS Low DNA is the premix PCR reagent developed by the use of Takara Bio's original technique for purification and DNA inactivation to suppress all detectable DNA in the reagent, derived from host cells and/or contaminated from environment, which can be template on PCR. In addition, this enables you to perform highly sensitive and specific PCR amplification quickly because of the use of modified Taq DNA polymerase having higher extension speed and specificity. This enzyme is designed for hot start PCR using the monoclonal antibody, thus it is possible to prevent non-specific amplification due to mis-priming and/or formation of primer dimer before PCR reaction. Since the antibody is denatured in the initial DNA-denaturation step, the additional activation step is not required. This enzyme has neither 5' -3' exonuclease nor 3' -5' exonuclease activities.

Storage: -20°C

Origin :

Escherichia coli carrying a plasmid that encodes the *Thermus aquaticus* DNA polymerase gene

Application :

- For DNA amplification of low copy template by PCR
- For microbial flora analysis

PCR products :

As most PCR products amplified with TaKaRa Taq HS Low DNA have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of them.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 20 µl) :

TaKaRa Taq HS Low DNA (2X)	10 µl
Template	< 100 ng
Primer 1	0.2 µM ^{*1} (final conc.)
Primer 2	0.2 µM ^{*1} (final conc.)
dH ₂ O for Low DNA	up to 20 µl

Note : Reaction mixture can be set up at room temperature. Be sure to keep all reagents on ice.

*1 Adjust the primer concentration to 0.6 µM, when you use high speed thermal cycler that ramp temperature speed is over 5°C.

Primer designing :

Each primer should be designed for T_m value more than 55°C. The T_m value should be calculated by below approximation calculation.

$$T_m = 4 \times (\text{the number of G or C bases}) + 2 \times (\text{the number of A or T bases}) + 35 - 2 \times (\text{total number of bases})$$

Note : In the case of calculating with Nearest Neighbor method, T_m value of designed primer should be more than 60°C.

PCR conditions :

2 step PCR

94°C ^{*2}	5 sec	} 30 - 35 cycles ^{*3}
65°C	20 sec/kb	

3 step PCR^{*4}

94°C ^{*2}	5 sec	} 30 - 35 cycles ^{*3}
55°C	1 sec	
68°C	20 sec/kb	

*2 Be sure to set the denaturation temperature at 94°C. Setting denaturation temperature at over 95°C might cause decrease of reactivity due to inactivation of enzyme.

*3 The setting of cycle number over 40 cycles is not recommended.

*4 In case of using primers which is designed for T_m value less than 55°C in approximation calculation (or less than 60°C in Nearest Neighbor method), select 3 step PCR reaction.

TaKaRa Taq is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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