

TECH NOTE

Comparison of the Guide-it Mutation Detection Kit with a CEL nuclease-based assay

Simple method to identify insertions or deletions in mammalian cells:

Amplify genomic regions directly from cells without the need for DNA purification, and detect mutations with the highly-efficient Guide-it Resolvase enzyme >>

Faster and more efficient than CEL nuclease-based assays:

The Guide-it mutation detection protocol is several hours shorter, more sensitive, and less prone to non-specific cleavage >>


Overview

Recently-developed genome editing tools such as zinc finger nucleases, transcription activator-like effector nucleases (TALENs), and the clustered regularly interspaced short palindromic repeats (CRISPR)/Cas9 system allow precise manipulation of virtually any gene. All of these editing techniques can be used to introduce double strand breaks at a target DNA sequence that are repaired by the error-prone non-homologous end joining (NHEJ) DNA repair pathway, resulting in introduction of insertion or deletion mutations (INDELS). Detecting these types of induced INDELS at target loci requires a simple and robust method.

A PCR-Based Method to Confirm the Presence of Mutations


Mutation detection is often based on PCR amplification of the region of interest and detection of mismatches in heteroduplexed DNA. With the [Guide-it Mutation Detection Kit](#), the target sequence is amplified directly from cells, without genomic DNA extraction/purification (Figure 1, step 1). Then, the PCR products are melted and rehybridized, forming mismatched targets that can be cleaved by the Guide-it Resolvase (Figure 1, steps 2 and 3).

1 Amplify genomic DNA from cells using Terra PCR Direct Polymerase




Region of interest is amplified (mutation is marked in blue)

2 Denature and reanneal



Guide-it Resolvase cleaves imperfectly matched DNA

3 Imperfectly matched DNA is cleaved by Guide-it Resolvase



4 Separate cleaved and uncleaved PCR fragments using agarose gel electrophoresis



Figure 1. Overview of the Guide-it Mutation Detection Kit method to confirm the presence of mutations in genomic DNA.

Comparison of the Guide-it Mutation Detection Kit with a CEL nuclease-based assay

The key component of the Guide-it Mutation Detection Kit is the Guide-it Resolvase, a mismatch-specific nuclease that recognizes heteroduplexed DNA. This enzyme is more efficient and more robust than other similar nucleases, such as Cel1.

To compare the Guide-it system and an assay based on CEL nuclease for detecting CRISPR/Cas9-introduced mutations in mammalian cells, 293T cells were transfected with plasmids encoding Cas9 and an sgRNA specific for the AAVS1 locus. Transfected cells harvested 48 hours post-transfection were mixed with untransfected cells at varying ratios (Figure 2, top). A DNA fragment containing the AAVS1 locus was generated by PCR using Terra Direct Polymerase, and the products were purified and cleaved with either Guide-it Resolvase (Guide-it Mutation Detection Kit) or the Cel1 enzyme (Company T). Mutations were easily discernible when using the Guide-it kit (Figure 2, bottom). In contrast, the CEL assay showed considerable smearing, making it difficult to determine cleavage efficiency and reducing the ability to detect lower levels of mutation (Figure 2, bottom).

	1	2	3	4	5	6
transfected	100	80	60	40	20	0
Non-transfected	0	20	40	60	80	100




Figure 2. Comparison of the Guide-it Mutation Detection kit and a CEL nuclease assay for detecting CRISPR/Cas9-introduced mutations in mammalian cells. Mutations were easily discernible when using the Guide-it kit (estimation of cleavage, 1: 59%, 2: 46%, 3: 28%, 4: 15%, 5: <10%, 6: 0%). In contrast, the CEL nuclease assay showed considerable smearing, making it difficult to determine cleavage efficiency and reducing the ability to detect low levels of mutation.

Summary

The Guide-it protocol, which amplifies genomic regions directly from cells without the need for DNA purification, reduced assay time from the 6 hours required by existing mismatch detection protocols to just 3.5 hours. In addition, this protocol provided increased mutation detection efficiency and demonstrated less sensitivity to buffers used in the PCR reaction than existing protocols.

In summary, compared to a CEL nuclease assay, the Guide-it mutation detection protocol is several hours shorter, more sensitive, and less prone to non-specific cleavage.

[Learn more about the Guide-it Mutation Detection Kit >>](#)

[Top of page](#)

[View web page >>](#)

http://www.clontech.com/US/Products/Cell_Biology_and_Epigenetics/Genome_Modification/CRISPR_Cas9/Guide-it_Resolvase_Versus_Cel1_Nuclease

This is a reprint from a page on our web site. All license, copyright, and trademark information pertaining to this content applies as stated in the original web content. This information can be found at www.clontech.com.