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v201003Da

High-Salt Solution for Precipitation (Plant)

Code No. 9193
Size : 50 ml

Shipping at RT
Stored at RT

Lot No.

* 2 years from date of receipt under proper storage conditions when unused. Once opened, it should be used earlier avoiding contamination.

Description : High-Salt Solution for Precipitation (Plant) is a reagent for removal of polysaccharides from RNA solution which extracted from plants containing a lot of polysaccharides and proteoglycans et al.. High purity of RNA sample is easily obtained only by addition of this product at isopropanol precipitation in purification steps of RNA extraction procedure.

Form :
1.2 M NaCl, 0.8 M Sodium Citrate

Usage :
Removal of polysaccharides from RNA solution.

Purity :
It is checking that this product is DNase and RNase free.

Protocol :

1. Add 0.5 volume of this solution to RNA sample*¹, and mix.
2. Add same volume of isopropanol as the solution added at 1), and mix.
3. Stand for 10 minutes at room temperature.
4. Centrifuge at 12,000 rpm for 10 minutes at 4 °C .
5. Discard the supernatant.
6. Add 1 ml of cold 75 % ethanol to the precipitate.
7. Re-centrifuge at 12,000 rpm for 5 minutes at 4 °C .
8. Discard the supernatant, and dry the precipitate briefly*².
9. Dissolve the precipitate in appropriate amount of RNase-free water or TE buffer.

* 1 : This product is also available for the extraction solution before isopropanol precipitation during RNA extraction procedure by AGPC method.

* 2 : Do not dry completely the precipitate. It may cause difficulty with dissolving RNA.


Notes :

1. It is not recommended to use for recovery of small molecular RNA, because of low recovery of that.
2. When RNA sample is small amount, it is recommended to use Dr. GenTLE™ Precipitation Carrier (Cat. #9094) at isopropanol precipitation for increasing the RNA recovery.

Application example :

Three samples of RNA solutions prepared from tobacco cultured cell (BY2), of which A₂₆₀/A₂₃₀ ratio are less than 1.5, are treated with this product for removal of polysaccharides according to the Protocol. After treatment, the A₂₆₀/A₂₃₀ ratio of the RNA samples became more than 2.0, showing that the polysaccharides in the RNA samples are removed.

Removal of polysaccharides

Related product :

Dr. GenTLE™ Precipitation Carrier (Cat. #9094)

Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.
Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.
If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at www.takara-bio.com .