# Brevibacillus Expression vector pNY326 DNA

Code No. HB111  
Size : 10  $\mu$ g (0.2 OD)

Shipping at  $-20^{\circ}\text{C}$   
Stored at  $-20^{\circ}\text{C}$

\* 2 years from date of receipt under proper storage conditions.

Lot No. PNY002

Concentration : 0.2  $\mu$ g /  $\mu$ l

Volume : 50  $\mu$ l

## Description :

pNY326 is an expression vector that can replicate in *B. choshinensis* cells. Because its size (3.4 kb) is smaller than that of pNCMO2 (5.2 kb) and the promoter activity of the vector is weaker than that of pNCMO2, this vector can be maintained stably, even if it expresses rather toxic proteins to the host cells. Even after repeated inoculation to scale up the culture, this vector still exhibits stable protein productivity. Thus, this vector is suited to culture on large scale. The vector can be maintained only in *B. choshinensis*. Because the construction of expression plasmid must be constructed in one step using *B. choshinensis* as the host, it is recommended to use the *Brevibacillus choshinensis* Electro-Cells ( $10^5$  transformants/ $\mu$ g DNA) as a transformation host.

**Form :** 10mM Tris-HCl, pH8.0  
1mM EDTA


**Base pairs :** 3,438 bp

## Usage :

This is a vector for protein expression using *Brevibacillus* Expression System.

**Cloning sites of pNY326 DNA :** (Please refer other side of this sheet)

## Map of pNY326 DNA :

## License agreement :

1. *B. choshinensis* SP3, pNY326, pNCMO2 and pNCMO2-BLA were developed and manufactured by Higeta Shoyu Co., Ltd. and sold by TAKARA BIO INC.
2. The use of this product is allowed for academic research only. You have to conclude a license agreement if you use this product or the products produced by using this product for commercial purpose. Contact us if you would like to ask any question about license etc. If any questions are received, either Higeta Shoyu Co., Ltd. or an agent for it will answer.  
TaKaRa Bio Inc.  
Tel +81-77-543-7247 Fax +81-77-543-9254
3. All rights of intellectual property related to this product are owned by Higeta Shoyu Co., Ltd. This product is covered by the following Japanese patents and corresponding international patent : Japanese patents, 2103774, 3696322, 3734593, 3753945, PCT/JP2004/016912 (international application based on PCT treaty), 3433807.
4. This product, the components or the derivatives of this product, or the products produced by using these items can not be transferred (distributed freely or sold) to third party. (When you purchase this product, you submit the license agreement to us.)

## Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.  
Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.  
If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at [www.takara-bio.com](http://www.takara-bio.com).

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