# T4 DNA Polymerase

**Code No. 2040A**      **Size:**            **100 U**  
**Conc.:**                **5 U/ $\mu$ l**

## Supplied Reagents:

**10X T4 DNA Polymerase Buffer**      **1 ml**  
**0.1% BSA**                                      **1 ml**

## Description:

T4 DNA polymerase, which requires template DNA and primer, selectively catalyzes the transfer of complementary dNTPs of the template to the 3'-OH terminus of the primer. This enzyme has 3' → 5' exonucleolytic activity specific for single-stranded DNA 10<sup>2</sup> - 10<sup>3</sup> - fold stronger than that of Klenow Fragment, but it does not have 5' → 3' exonuclease activity.

## Storage Buffer:

200 mM	Potassium Phosphate, pH 6.5
1 mM	DTT
50%	Glycerol

**Storage:**            -20°C

## Source:

*Escherichia coli* carrying the plasmid containing phage T4 DNA polymerase gene.

## Unit definition:

One unit is the amount of the enzyme catalyzing the incorporation of 10 nmol of total nucleotides into acid-insoluble products in 30 min at 37°C and pH 8.8, with heat-denatured calf thymus DNA as the template-primer.

## Reaction mixture for unit definition:

67 mM	Tris-HCl, pH 8.8
6.7 mM	MgCl <sub>2</sub>
16.6 mM	(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>
1 mM	DTT
6.7 $\mu$ M	EDTA
0.0167%	bovine serum albumin
0.2 mg/ml	heat-denatured calf thymus DNA
330 $\mu$ M each	dCTP, dATP, and dGTP
330 $\mu$ M	[ <sup>3</sup> H] dTTP

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. Making the termini of DNA fragments blunt-ended.<sup>4)</sup>
2. Labeling from the 3'-termini of DNA fragments by replacement synthesis with the use of its strong 3' → 5' exonuclease activity.<sup>3)</sup>
3. Analysis of the starting point of transcription of mRNA by the primer extension method.<sup>5)</sup>

## Note:

The optimum pH is 8 - 9; at pH 7.5 and 9.7, the activity decreases by 50%. The enzyme requires Mg<sup>2+</sup> and SH reagents for maximum activity.<sup>1)</sup> The enzyme is inhibited when the total ionic strength of the reaction mixture exceeds 100 mM.<sup>2)</sup> This enzyme tends to be affected by the steric structure of the template DNA, such as in having its polymerase activity activated by the T4 protein coded by gene 32, at which time the 3' → 5' exonuclease activity is completely inhibited.

## Composition of Supplied Reagents (Store at -20°C):

1. 10X T4 DNA Polymerase Buffer  
330 mM Tris-acetate, pH 7.9  
660 mM Potassium acetate  
100 mM Magnesium acetate  
5 mM DTT

2. 0.1% BSA

The supplied buffer has the different composition from that used for unit definition. The supplied buffer has the general composition to be applicable to blunting of the end of double-stranded DNA.

White precipitate may be generated when adding 0.1% BSA directly to the supplied buffer. Prepare the reaction mixture (final conc. 0.01% BSA) by adding the reagents in the order given below.

Sterile purified water → supplied buffer → 0.1% BSA → substrate DNA

## Application example:

### Blunting the termini of DNA fragments

1. Prepare the reaction mixture by combining the followings in a microcentrifuge tube to the total volume of 9  $\mu$ l.

Insert DNA with protruding end	> 0.1 pmol
10X T4 DNA Polymerase Buffer	1 $\mu$ l
0.1% BSA	1 $\mu$ l
1.7 mM dNTP Mixture	1 $\mu$ l
Sterile purified water	up to 9 $\mu$ l
2. To avoid annealing of the termini of DNA, heat at 70°C for 5 min and place at the 37°C incubation bath.
3. Add 1  $\mu$ l of T4 DNA Polymerase (1 U/ $\mu$ l : dilute with 1X supplied buffer etc.) and mix gently by pipetting. Avoid vigorous mixing with Vortex mixer.
4. Incubate at 37°C for 5 min.
5. Mix vigorously with a Vortex mixer to inactivate the enzyme. (Almost no enzyme activity remains after Vortex mixing. Place on ice to ensure the inactivation of T4 DNA polymerase and immediately perform the next ligation. If you wish to stop the reaction before ligation, the mixture may be stored at -20°C following phenol/chloroform extraction and ethanol precipitation.)

## References:

- 1) Lehman I R. *Methods in Enzymology*. (1974) **29**: 46-53.
- 2) Goulian M, Lucas Z J, and Kornberg A. *J Biol Chem*. (1968) **243**: 627-638.
- 3) Deen K C, Landers T A, and Berninger M. *Anal Biochem*. (1983) **135**: 456-465.
- 4) Wartell R M, and Reznikoff W S. *Gene*. (1980) **9**: 307-319.
- 5) Hu M C T and Davidson N. *Gene*. (1986) **42**: 21-29.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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