Ribonuclease H (RNase H)

Code No. 2150A **Size:** **1,000 U**
Conc.: **60 U/μl**

Description :

RNase H catalyzes the specific degradation of RNA in DNA-RNA hybrids. Its MW is 21,000 and its optimum pH is about 8.0. It is activated by Mg²⁺ or Mn²⁺.

Storage Buffer :

25 mM	Tris-HCl, pH 7.5
30 mM	NaCl
0.5 mM	EDTA
1 mM	DTT
50%	Glycerol

Storage: -20°C

Source :

Escherichia coli HB101 containing *rnh* plasmid (pKH11) and regulator plasmid (pNT203)¹⁻⁴

Unit definition :

One unit is the amount of the enzyme that produces 1 nmol of acid-soluble products in 20 minutes at 30°C and pH 7.7, with [³H] poly (rA) · poly (dT) as the substrate.

Reaction mixture for unit definition :

40 mM	Tris-HCl, pH 7.7
4 mM	MgCl ₂
1 mM	DTT
4 %	Glycerol
0.003 %	bovine serum albumin
24 μM	[³ H] poly (rA) · poly (dT)

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Note :

This enzyme is obtained from an *E. coli* overproducing strain and is pure, containing very little contamination. Its concentration is high, so problems that might arise if an excess of the enzyme is used should not occur.

Applications :

1. cDNA cloning by the Okayama-Berg method⁵
2. Detection of DNA-RNA hybrids⁶
3. Removal of the poly(A) tail from mRNA in the presence of oligo(dT)⁷

References :

- 1) Berkower I, Leis J, and Hurwitz J. *J Biol Chem.* (1973) **248**: 5914-5921.
- 2) Horiuchi T, Maki H, Maruyama M, and Sekiguchi M. *Proc Natl Acad Sci USA.* (1981) **78**: 3770-3774.
- 3) Maki H, Horiuchi T, and Sekiguchi M. *Proc Natl Acad Sci USA.* (1983) **80**: 7137-7141.
- 4) Shigesada K, Tsurushita N, Matsumoto W, and Imai M. *Gene.* (1984) **29**: 199-209.
- 5) Okayama H and Berg P. *Mol Cell Biol.* (1982) **2**: 161-170.
- 6) Keller W and Crouch R. *Proc Natl Acad Sci USA.* (1972) **69**: 3360-3364.
- 7) Vournakis J N, Efstratiadis A, and Kafatos F C. *Proc Natl Acad Sci USA.* (1975) **72**: 2959-2963.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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