# Recombinant DNase I (RNase-free)

**Code No. 2270A**      **Size:** 1,000 U  
**Conc.:** 5 U/ $\mu$ l

**Supplied Reagent:**  
**10X DNase I Buffer**      **1 ml**

## Description:

Recombinant DNase I (RNase-free) is an endodeoxyribonuclease that catalyzes to the same extent the degradation of both single- and double-stranded DNA randomly, and produces 5'-P terminal oligonucleotides. As protease activity and ribonuclease activity have been eliminated, this enzyme is stable around its optimum neutral pH range. So, this enzyme is suitable for RNA preparation at neutral pH.

## Storage Buffer:

20 mM	Tris-HCl, pH 7.5 (at 4 °C)
50 mM	NaCl
0.1 mM	CaCl <sub>2</sub>
50 %	Glycerol

**Storage:** -20°C

**Source:** Recombinant enzyme derived from non-animal host.

## Unit definition:

One unit is the amount of the enzyme that increases the absorbance at 260 nm by 0.001 per minute at 25°C, pH 5.0, with calf thymus DNA as the substrate (Kunitz unit).<sup>1)</sup>

## Reaction mixture for unit definition:

100 mM	Sodium acetate, pH 5.0 (at 25 °C)
5 mM	MgSO <sub>4</sub>
40 $\mu$ g/ml	calf thymus DNA

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. For digestion of template DNA without buffer exchange, after *in vitro* transcription with T7 or SP6 RNA Polymerase.<sup>5)</sup>
2. For nick translation together with DNA Polymerase I (Cat.#2130).<sup>2)</sup>
3. Making a DNA library for shotgun sequencing in the presence of Mn<sup>2+</sup>.<sup>3)</sup>
4. For foot-printing analysis of the interaction between DNA and protein.<sup>4)</sup>
5. For digestion of genomic DNA before RT-PCR.

## Note:

As protease activity has been eliminated, this enzyme is stable around its optimum neutral pH range. The enzyme requires a divalent metal ion for its activity. It randomly produces nicks in double-stranded DNA in the presence of Mg<sup>2+</sup>, but in the presence of Mn<sup>2+</sup>, both strands of double stranded DNA are cleaved into fragments. The enzyme loses its activity reversibly with EDTA, and irreversibly by heat treatment at 75°C for 10 min.

## Composition of Supplied Reagents (Store at -20 °C):

10X DNase I Buffer	
400 mM	Tris-HCl, pH 7.5
80 mM	MgCl <sub>2</sub>
50 mM	DTT

## Application example:

### Digestion of genomic DNA in a sample for RT-PCR

1. Prepare the following reaction mixture.

total RNA	20 - 50 $\mu$ g
10X DNase I Buffer	5 $\mu$ l
Recombinant DNase I (RNase-free)	2 $\mu$ l (10 U)
Ribonuclease Inhibitor	20 U
DEPC-treated water	up to 50 $\mu$ l
2. Incubate for 20 - 30 min at 37°C.
3. Perform one of the following procedures to inactivate DNase I.
  - A. Heat treatment
    - (1) Add 2.5  $\mu$ l of 0.5 M EDTA, and incubate at 80°C for 2 min.
    - (2) Increase reaction volume to 100  $\mu$ l with DEPC treated water.
  - B. Phenol/Chloroform extraction
    - (1) Mix 50  $\mu$ l of DEPC treated water and 100  $\mu$ l of phenol/chloroform/isoamyl alcohol (25 : 24 : 1) together.
    - (2) Centrifuge at 12,000 rpm for 5 min at room temperature, then transfer the upper layer to a new tube.
    - (3) Add equal amount of chloroform/isoamyl alcohol (24 : 1) and mix.
    - (4) Centrifuge at 12,000 rpm for 5 min at room temperature, then transfer upper layer to new tube.
4. Add 10  $\mu$ l of 3 M sodium acetate and 250  $\mu$ l of chilled ethanol, and then mix. Keep it for 20 min at -80°C.
5. Centrifuge at 12,000 rpm for 10 min at 4 °C. Remove the supernatant.
6. Wash the precipitate with chilled 70 % ethanol. Centrifuge at 12,000 rpm for 5 min at 4 °C and remove the supernatant.
7. Dry the precipitate.
8. Dissolve the precipitate in a suitable amount of DEPC-treated water. Confirm the genomic DNA is removed completely by electrophoresis and measure the RNA concentration. When the genomic DNA is not removed completely, increase the amount of enzyme or extend reaction time.

## References:

- 1) Kunitz M. *J Gen Physiol.* (1950) **33**: 349-377.
- 2) Rigby P W J, Dieckmann M, Rhodes C, and Berg P. *J Mol Biol.* (1977) **113**: 237-251.
- 3) Anderson S. *Nucleic Acids Res.* (1981) **9**: 3015-3027.
- 4) Galas D J and Schmitz A. *Nucleic Acids Res.* (1978) **5**: 3157-3170.
- 5) Green M R, Maniatis T, and Melton D A. *Cell.* (1983) **32**: 681-694.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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● 起源      Recombinant enzyme derived from non-animal host.

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