# Recombinant RNase Inhibitor ver.2.0

**Code No. 2315A**      **Size:**      **5,000 U**  
**Conc.:**      **40 U/μl**

## Description:

Recombinant RNase Inhibitor ver.2.0 is an updated product of the recombinant protein of porcine liver RNase inhibitor. While this product maintains characteristics similar to RNase inhibitors from porcine liver and human placenta,<sup>1,2)</sup> the stability of this protein has been increased by introducing a mutation in one of the cysteine residues, which is responsible for the denaturation due to its easy oxidation.<sup>3,4)</sup> It forms a 1 : 1 complex with RNase A to inhibit RNase activity.<sup>5)</sup> This reaction is reversible, and the inhibitor can be irreversibly inactivated to restore ribonuclease activity by dissociating the complex with urea or sulfhydryl reagent. Moreover, unlike other non-protein competitive inhibitors (e.g., nucleotides and inorganic phosphates), it can easily be removed from the reaction system by phenol extraction. In addition, this product does not inhibit RNase H activity. This product can be directly added in various reaction mixtures such as *in vitro* transcription and RT-PCR, where the integrity of RNA is indispensable.

**Storage Buffer:**      20 mM    HEPES-KOH, pH 7.5  
                                 50 mM    KCl  
                                 5 mM    DTT  
                                 50%    Glycerol

**Storage:**      -20°C

## Source:

*E. coli* carrying the plasmid containing the gene for ribonuclease inhibitor from porcine liver

## Properties:

Molecular mass : approx. 52 kDa  
Optimum pH : Maximum at pH 7 - 8, though the inhibition is observed in a broad pH range

## Unit Definition:

One unit is defined as the amount of protein that inhibits 50% of RNase A activity, when 2',3'-cyclic CMP is used as the substrate for 5 ng of RNase A.

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Precaution for Use:

1. Don't mix the protein vigorously.
2. Include some DTT in the reaction to secure the activity, as far as its concentration doesn't interfere with the reaction (e.g., 1 mM).

## Applications:

1. *In vitro* transcription/translation (1 U/μl reaction)<sup>7,8)</sup>
2. *In vitro* transcription/translation with cell-free extract (20 U/μl reaction)<sup>7)</sup>
3. RT-PCR (0.5 U/μl reaction)
4. cDNA synthesis (0.5 U/μl reaction)<sup>9)</sup>
5. Polysome isolation (1 U/μl reaction)<sup>8)</sup>

\* The parentheses show the examples of RNase inhibitor concentration in each reaction mixture.

## References:

- 1) Burton L E and Fucci N P. *Int J Pept Protein Res.* (1982) **19**: 372-379.
- 2) Blackburn P, Wilson G, and Moore S. *J Biol Chem.* (1977) **252**: 5904-5910.
- 3) Kim B-M, Schultz L W, and Raines R T. *Protein Sci.* (1999) **8**: 430-434.
- 4) Dickson K A, Haigis M C, and Raines R T. *Prog Nucleic Acid Res Mol Biol.* (2005) **80**: 349-374.
- 5) Turner P M, Lerea K M, and Kull F J. *Biochem Biophys Res Commun.* (1983) **114**: 1154-1160.
- 6) Blackburn P. *J Biol Chem.* (1979) **254**: 12484-12487.
- 7) Eichler D C, Tatar T F, and Lasater L S. *Biochem Biophys Res Commun.* (1981) **101**: 396-403.
- 8) Scheele G and Blackburn P. *Proc Natl Acad Sci USA.* (1979) **76**: 4898-4902.
- 9) de Martynoff G, Pays E, and Vassart G. *Biochem Biophys Res Commun.* (1980) **93**: 645-653.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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1. *In vitro* transcription/translation (1 U/ $\mu$ l reaction)<sup>7,8)</sup>
2. *In vitro* transcription/translation with cell-free extract (20 U/ $\mu$ l reaction)<sup>7)</sup>
3. RT-PCR (0.5 U/ $\mu$ l reaction)
4. cDNA synthesis (0.5 U/ $\mu$ l reaction)<sup>9)</sup>
5. Polysome isolation (1 U/ $\mu$ l reaction)<sup>8)</sup>

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