Recombinant RNase Inhibitor, HQ

Code No. 2316A **Size:** **40,000 U**
Conc.: **40 U/ μ l**

Description:

Recombinant RNase Inhibitor, HQ (high quality) is intended to be used for the preparation of active pharmaceutical ingredients in non-clinical trials and for the process development of drug manufacturing under GMP regulation. It can also be used for basic research of RNA therapeutics. The product is free from not only any human- or animal-derived materials but also β -lactam compounds in the final formulation. This product is a variant of the ribonuclease inhibitor derived from porcine liver and has an equivalent performance to Recombinant RNase Inhibitor ver.2.0 (Cat. #2315A/B). While this product maintains characteristics similar to RNase inhibitors from porcine liver and human placenta,^{1,2} the stability of this protein has been increased by introducing a mutation in one of the cysteine residues, which is responsible for the denaturation due to its easy oxidation.^{3,4} This product can be directly added in various reaction mixtures such as *in vitro* transcription and capping reaction, where the integrity of RNA is indispensable. Therefore, the product is suitable for research and development in the field of RNA therapeutics such as vaccines.

Product Quality:

1. This product does not contain any human- or animal-derived materials in the final formulation.
2. This product does not contain any β -lactam compounds in the final formulation.

Storage Buffer:

20 mM	HEPES-KOH, pH 7.5
50 mM	KCl
5 mM	DTT
50%	Glycerol

Storage: -20°C

Source:

Escherichia coli carrying a plasmid containing the gene for a variant of ribonuclease inhibitor from porcine liver

Properties:

Molecular mass : approx. 52 kDa
Optimum pH : Maximum at pH 7 - 8, though the inhibition is observed in a broad pH range

Unit Definition:

One unit is defined as the amount of protein that inhibits 50% of RNase A activity, when 2',3'-cyclic CMP is used as the substrate for 5 ng of RNase A.⁵

Reaction Mixture for Unit Definition:

100 mM	Tris-acetate pH 6.5 (26°C)
1 mM	EDTA
1 mM	2',3'-Cyclic CMP
1 μ g/ml	RNase A

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. *In vitro* transcription/translation (1 U/ μ l reaction)^{6,7}
2. *In vitro* transcription/translation with cell-free extract (20 U/ μ l reaction)⁶
3. Capping reaction (1 U/ μ l reaction)

* The parentheses show the examples of RNase inhibitor concentration in each reaction mixture.

Precautions for Use:

1. Don't mix the protein vigorously.
2. Include DTT in the reaction to maintain the activity, as far as its concentration doesn't interfere with the reaction (e.g., 1 mM).

References:

- 1) Burton L E and Fucci N P. *Int J Pept Protein Res.* (1982) **19**: 372-379.
- 2) Blackburn P, Wilson G, and Moore S. *J Biol Chem.* (1977) **252**: 5904-5910.
- 3) Kim B-M, Schultz L W, and Raines R T. *Protein Sci.* (1999) **8**: 430-434.
- 4) Dickson K A, Haigis M C, and Raines R T. *Prog Nucleic Acid Res Mol Biol.* (2005) **80**: 349-374.
- 5) Blackburn P. *J Biol Chem.* (1979) **254**: 12484-12487.
- 6) Eichler D C, Tatar T F, and Lasater L S. *Biochem Biophys Res Commun.* (1981) **101**: 396-403.
- 7) Scheele G and Blackburn P. *Proc Natl Acad Sci USA.* (1979) **76**: 4898-4902.

Related Products:

Note: For more information on the "mRNA synthesis learning center", please see our website (<https://www.takarabio.com/>).

[HQ grade]

PrimeCap™ T7 RNA Polymerase (low dsRNA), HQ (Cat. #2561A)
T7 RNA Polymerase, HQ (Cat. #2542A)
Pyrophosphatase (inorganic), HQ (Cat. #2451A)
Faustovirus Capping Enzyme (S17), HQ (Cat. #2481A)
mRNA Cap 2'-O-Methyltransferase, HQ (Cat. #2471A) etc.

[RUO grade]

Recombinant RNase Inhibitor ver.2.0 (Cat. #2315A/B)
PrimeCap™ T7 RNA Polymerase (low dsRNA) (Cat. #2560A)
T7 RNA Polymerase ver.2.0 (Cat. #2541A) etc.

PrimeCap is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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