Vaccinia Capping Enzyme

Code No. 2460A **Size: 500 U**
Conc.: 10 U/μl

Supplied Reagents:

10X Capping Buffer **100 μl**
S-adenosylmethionine (SAM) (32 mM) **100 μl**
GTP (10 mM) **50 μl**

Description:

Vaccinia Capping Enzyme is a capping enzyme originating from vaccinia virus and adds a 7-methylguanylate structure (Cap 0) to the 5' end of RNA (5'-triphosphate RNA). In eukaryotes, cap structures are important for stabilization, nuclear export, and translation of mRNAs. This enzyme consists of two subunits (D1 and D12) and has three enzymatic activities (RNA triphosphatase, guanylyltransferase, and guanine methyltransferase). The use of this enzyme together with the reagents included in this product (10X Capping Buffer, S-adenosylmethionine (SAM), and GTP) allows to produce an intact Cap-0 RNA efficiently *in vitro*. Moreover, the use of this enzyme together with mRNA Cap 2'-O-Methyltransferase (Cat. #2470A/B) allows to produce an intact Cap-1 RNA that have a methyl group at the 2'-O position of the first nucleotide of Cap-0 RNA by a single-step reaction.

Storage: -20°C

Source:

Escherichia coli carrying plasmids containing the genes for vaccinia virus capping enzyme

Properties:

Molecular mass: approx. 98 kDa and 33 kDa (heterodimer)

Unit Definition:

One unit is defined as the amount of enzyme that incorporates 10 pmol of [α -32P] GTP into an RNA transcript in 1 hour at 37°C.

Reaction Mixture for Unit Definition:

1X	Capping Buffer
1 μM	GTP
10 μCi/ml	[α -32P] GTP
0.67 mM	SAM
2.5 μg/50 μl	100 nt RNA

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. Preparation of Cap-0 mRNA
2. Preparation of 5' end-labeled mRNA

Precautions for Use:

1. Don't mix the enzyme vigorously.
2. RNA is less produced or fragmented when reagents, tubes, or pipette tips are contaminated with RNase. Precautions should be taken during experiments to avoid RNase contamination, such as wearing disposable gloves and using tubes and tips exclusively for RNA experiments.

Protocols:

1. Preparation of Cap-0 RNA from uncapped RNA (5'-triphosphate RNA)

Denatured RNA transcript (10 μg)*1,4	15 μl
10X Capping Buffer	2 μl
GTP (10 mM)	1 μl
SAM (2 mM, dilute 32 mM stock to 2 mM)*2	1 μl
Vaccinia Capping Enzyme (10 U/μl)	1 μl
<hr/>	
Total	20 μl*3

Incubate at 37°C for 30 minutes. *4

2. Preparation of Cap-1 RNA from uncapped RNA (5'-triphosphate RNA)

Denatured RNA transcript (10 μg)*1,4	14 μl
10X Capping Buffer	2 μl
GTP (10 mM)	1 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	1 μl
Vaccinia Capping Enzyme (10 U/μl)	1 μl
mRNA Cap 2'-O-Methyltransferase (50 U/μl)	1 μl
<hr/>	
Total	20 μl*3

Incubate at 37°C for 1 hour. *4

- *1 In order to remove secondary structure on the 5' end of the transcript, heat-denature the RNA (5' triphosphate RNA) prior to the reaction.
 - (1) Take 10 μg of the RNA and adjust the volume to 15 μl (Cap-0 RNA prep) or 14 μl (Cap-1 RNA prep) by RNase-free water.
 - (2) Heat-denature at 65°C for 5 - 10 min, immediately followed by an incubation on ice for 5 min.
- *2 SAM is unstable. Dilute necessary amount of 32 mM stock solution by RNase-free water prior to the reaction, and keep it on ice until use.
- *3 It can be scaled up depending on the amount required.
- *4 If the RNA length is less than 200 nt, extend the incubation time to enhance the capping efficiency (and methylation efficiency) (Cap-0 RNA prep: 30 min → 1 hr, Cap-1 RNA prep: 1 hr → 2 hrs).

Note: RNase inhibitor can be used to enhance the stability of RNA during the reaction. Add at a final concentration of 1 U/μl if necessary.

References:

- 1) Shuman S. *J Biol Chem.* (1990) **265**: 11960-11966.
- 2) Furuichi Y, LaFiandra A, and Shatkin A J. *Nature.* (1977) **266**: 235-239.
- 3) Lewis J D and Izaurralde E. *Eur J Biochem.* (1997) **247**: 461-469.
- 4) Ramanathan A, Bobb G B, and Chan S H. *Nucleic Acids Res.* (2016) **44**: 7511-7526.

Related Products:

Takara IVTpro™ mRNA Synthesis System (Cat. #6141)
Cloning Kit for mRNA Template (Cat. #6143)
Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)
T7 RNA Polymerase ver.2.0 (Cat. #2541A)
mRNA Cap 2'-O-Methyltransferase (Cat. #2470A/B)
Recombinant RNase Inhibitor (Cat. #2313A/B)
RNase-free Water (Cat. #9012)
NucleoSpin RNA Clean-up (Cat. #740948.10/.50/.250)*

* Not available in all geographic locations. Check for availability in your area.

IVTpro is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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