NucleoGone™ Endonuclease

Code No. 2490A **Size: 25,000 U**
Conc.: 100 U/μl

Description:

NucleoGone Endonuclease is an enzyme that efficiently degrades DNA and RNA (single-stranded, double-stranded, linear, and circular forms) and works in a wide range of conditions. It is useful for degrading and removing of existing nucleic acids in processes such as protein purification and viral vector purification from host cells.

Storage: -20°C

Source:

Escherichia coli carrying a plasmid that encodes the gene for NucleoGone Endonuclease

Property:

Molecular mass: approx. 27 kDa

Unit Definition:

One unit is defined as the amount of enzyme that changes the absorbance at 260 nm by 1.0 in 30 minutes at 37°C, using activated salmon sperm DNA as the substrate.

Reaction Mixture for Unit Definition:

50 mM	Tris-HCl, pH 8.0
1 mM	MgCl ₂
0.01%	BSA
0.5 mg/ml	activated salmon sperm DNA

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. Removal of genomic DNA derived from host cells during virus purification, such as AAV.
2. Reduction of viscosity of cell lysates during protein purification.

Optimal Conditions, Inhibition, and Inactivation:

1. NucleoGone Endonuclease retains its activity under the following conditions.

Condition	Optimal	Effective
Temperature	37 - 42°C	4 - 42°C
pH	8 - 10	7 - 10
Na ⁺ , K ⁺	0 - 100 mM	0 - 300 mM
Mg ²⁺	≥ 1 mM	≥ 0.01 mM
DTT	0 - 50 mM	≥ 0 mM
2-mercaptoethanol	0 - 200 mM	≥ 0 mM
Sodium deoxycholate	0 - 0.1%	0 - 1%
Triton X-100	0 - 1%	0 - 1%

2. The activity is completely inhibited with ≥ 2 mM EDTA.
3. To completely inactivate the enzyme, add NaOH to 100 mM and heat at 70°C for 30 minutes.

NucleoGone is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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