# Reverse Transcriptase XL (AMV)

Code No. 2620A      Size:      500 U

## Supplied Reagent: 10X/5X Reaction Buffer

### Storage:

In order to conserve optimal activity, it is strongly suggested that the enzyme be aliquoted and stored at -80°C, although storage at -20°C is adequate for short periods of time. Repeated freezing and thawing results in loss of enzyme activity. Do not store this enzyme in a frost-free freezer.

### Description:

Reverse Transcriptase XL (AMV) is a RTase that is isolated from Avian myeloblastosis virus. It consists of the  $\alpha$   $\beta$  holoenzyme of molecular weight 157,000 daltons. This enzyme is highly purified, free of nucleases, and is qualified for use in cDNA synthesis using RNA templates approaching 10 kb in length. This enzyme has an associated RNase H activity independent of cDNA synthetic activity.

### Storage Buffer:

200 mM	Potassium Phosphate, pH 7.2
2 mM	DTT
0.2 %	Triton X-100
50 %	Glycerol (v/v)

**Source:** Avian myeloblastosis virus

### Unit definition:

One unit is defined as the amount of RTase XL (AMV) required to catalyze the incorporation of one nmol of dTTP into an acid-insoluble product in 10 minutes at 37°C using poly (rA) · (dT)<sub>12-18</sub> as template primer.

### Assay Conditions for unit definition:

50 mM	Tris-HCl, pH 8.3
40 mM	KCl
6 mM	MgCl <sub>2</sub>
4 mM	DTT
0.4 mM	poly (rA) · (dT) <sub>12-18</sub>
0.5 mM	[ <sup>3</sup> H]dTTP
2 - 4 U	RTase XL (AMV)

### Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

### Note :

The optimum temperature range appears to be 42 - 58°C. When performing cDNA synthesis at temperature higher than 50°C, a short pre-incubation at 42°C for 4 minutes is recommended.

### Composition of Supplied Reagent:

10X/5X Reaction Buffer	
250 mM	Tris-HCl, pH 8.3
500 mM	KCl
20 mM	DTT
50 mM	MgCl <sub>2</sub>

\* It is recommended to use this buffer at 10 - fold dilution for usual cDNA synthesis.

### References:

- 1) Hellman GM, Shahabuddin M, Shaw JG, and Rhoads RE. *Virology*. (1983) **128**: 210-220.
- 2) Houts GE, Miyagi M, Ellis C, Beard D, and Beard J W. *J Virol*. (1979) **29**: 517-522.
- 3) Shimomaye E and Salvato M. *Gene Anal Tech*. (1989) **6**: 25-28.
- 4) Fuchs B, Zhang K, Rock MG, Bolander ME, and Sarkar G. *Mol Biotechnol*. (1999) **12**: 237-40.
- 5) Maniatis T, *et al*. Cold Spring Harbor Laboratory, *Molecular Cloning*. (1982) 217-229.

This product is manufactured by Life Sciences Advanced Technologies, Inc.

### Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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