# PrimeScript™ Reverse Transcriptase

**Code No. 2680A**      **Size:**      **10,000 U**  
**Conc.:**      **200 U/μl**

**Supplied Reagents:**  
**5X PrimeScript Buffer**      **500 μl**

## Description:

PrimeScript Reverse Transcriptase is a modified M-MLV (Moloney Murine Leukemia Virus) RTase. This enzyme has extremely high extension capability and can synthesize long 1st-strand cDNA efficiently. Even for difficult templates, including RNAs with complex secondary structure, it is possible to synthesize 1st-strand cDNA at the normal reverse transcription temperature (42°C) with this enzyme. It is not necessary to perform the RT reaction at higher temperatures, a condition that may cause RNA degradation. This enzyme is suitable for preparation of long cDNAs, construction of cDNA libraries that include full-length cDNA, etc.

## Storage buffer:

20 mM	Tris-HCl, pH 7.8
100 mM	NaCl
1 mM	EDTA
1 mM	DTT
50%	Glycerol (v/v)

**Storage:**      -20°C

## Source:

*E. coli* carrying the plasmid that encodes reverse transcriptase gene.

## Unit definition:

One unit is the amount of the enzyme that incorporates 1 nmol of [<sup>3</sup>H] dTTP in 10 minutes at 37°C, with poly (rA) · oligo (dT)<sub>12-18</sub> as the primer-template.

## Reaction mixture for unit definition:

50 mM	Tris-HCl, pH 8.3
75 mM	KCl
8 mM	MgCl <sub>2</sub>
10 mM	DTT
20 μg/ml	(rA) <sub>n</sub> · (dT) <sub>12-18</sub>
0.5 mM	[ <sup>3</sup> H] dTTP
0.1%	NP-40

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. First-strand cDNA synthesis.
2. Preparation of cDNA probes.
3. RT-PCR.

## Composition of Supplied Reagent:

5X PrimeScript Buffer (for cDNA synthesis)  
250 mM Tris-HCl, pH 8.3  
375 mM KCl  
15 mM MgCl<sub>2</sub>

## Standard protocol for 1st-strand cDNA synthesis:

1. Prepare the following mixture in a microtube.

Oligo dT primer	50 pmol
(or random primer (6 mers))	50 pmol)
(or gene specific primer	2 pmol)
dNTP Mixture (10 mM each)	1 μl
Template RNA	total RNA ≤ 5 μg, mRNA ≤ 1 μg
RNase free dH <sub>2</sub> O	up to 10 μl

2. Heat at 65°C for 5 min and cool immediately on ice.
3. Prepare the reaction mixture by combining the following in a total volume of 20 μl.

Template RNA/Primer mixture	10 μl
5X PrimeScript Buffer	4 μl
RNase Inhibitor	20 U
PrimeScript Reverse Transcriptase	100 - 200 U *1
RNase free dH <sub>2</sub> O	up to 20 μl

\* 1 100 U is recommended for RT-PCR and cDNA cloning;  
200 U is recommended for quantitative analysis, such as qPCR.

4. Mix gently.
5. Perform the reaction under the following condition.  
30°C      10 min \*1  
42 (- 50)°C \*2      30 - 60 min \*3  
\* 1 This step is required for random primer.  
\* 2 It is generally recommended to perform the RT reaction at 42°C. However, for RT-PCR, if the reverse primer for PCR is also used as a RT primer, non-specific products may be amplified due to mispriming. In such a case, perform the RT reaction at 50°C for 30 min.  
\* 3 In most cases, 30 min is sufficient. Increase the incubation time to 60 min when the target is very long.
6. Heat at 70°C for 15 min and cool on ice.

The obtained cDNA can be used for 2nd-strand cDNA synthesis or as a template for PCR.

PrimeScript is a trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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