PrimeScript™ II Reverse Transcriptase

Code No. 2690A Size: 10,000 U
Conc.: 200 U/μl

Supplied Reagent:
5X PrimeScript II Buffer 500 μl

Description :

PrimeScript II Reverse Transcriptase (RTase) yields high-quality cDNA with the least amount of background at standard reaction temperature (42°C) that doesn't cause RNA degradation because of its prominent primer use efficiency and addition of accessory protein, which resolves RNA secondary structure to inhibit cDNA synthesis and prevents RT-primer dependent mis-priming. The accessory protein completely inhibits non-specific extension reaction caused by reverse transcriptase on ice, thus the inhibition of long-strand cDNA synthesis doesn't occur even if it takes a while to start reverse transcription reaction after the reaction mixture was prepared. The synthesized cDNA can be used for many experiments such as 2nd-strand synthesis, hybridization, endpoint and quantitative PCR, and also this enzyme is particularly suitable for applications that require high-quality cDNA such as construction of full-length cDNA library and so on.

Storage buffer :

20 mM	Tris-HCl, pH 7.8
100 mM	NaCl
1 mM	EDTA
1 mM	DTT
1%	Tween 20 (v/v)
50%	Glycerol (v/v)

Storage : -20°C

Source :

Purified from an *E. coli* strain expressing a recombinant enzyme.

Unit definition :

One unit is the amount of the enzyme that incorporates 1 nmol of [³H] dTTP in 10 minutes at 37°C, with poly (rA) · oligo (dT)₁₂₋₁₈ as the primer-template.

Reaction mixture for unit definition :

50 mM	Tris-HCl, pH 8.3
75 mM	KCl
8 mM	MgCl ₂
10 mM	DTT
20 μg/ml	(rA) _n · (dT) ₁₂₋₁₈
0.5 mM	[³ H] dTTP
0.1 %	NP-40

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

First-strand cDNA synthesis

- Preparation of cDNA probe
- RT-PCR

Composition of Supplied Reagent :

5X PrimeScript II Buffer (for cDNA synthesis)	
250 mM	Tris-HCl, pH 8.3
375 mM	KCl
15 mM	MgCl ₂

Standard protocol for 1st-strand cDNA synthesis :

1. Prepare the following mixture in a microtube.

Oligo dT primer	50 pmol
(or Random primer (6 mers)	50 pmol (20 - 100 pmol ^{*1}))
(or Gene specific primer	2 pmol)
dNTP Mixture (10 mM each)	1 μl
Template RNA	total RNA ≤ 5 μg ^{*2} , mRNA ≤ 1 μg
RNase free dH ₂ O	up to 10 μl

2. Heat at 65°C for 5 min and cool immediately on ice.

3. Prepare the reaction mixture by combining the following to a total volume of 20 μl.

Template RNA/Primer mixture	10 μl
5X PrimeScript II Buffer	4 μl
RNase Inhibitor	20 U
PrimeScript II RTase	200 U
RNase free dH ₂ O	up to 20 μl

4. Mix gently.

5. Perform the reaction under the following condition.

30°C	10 min (Only when Random 6 mers is used.)
42 (- 50)°C ^{*3}	30 - 60 min

6. Heat at 70°C for 15 min and cool on ice.

*1 For the use of random primer, use 20 - 50 pmol in case of cDNA synthesis for strand length more than 2 kb, on the other hand, use 50 - 100 pmol in case of cDNA synthesis for quantitative PCR.

*2 The amount of total RNA should be less than 1 μg in case of cDNA synthesis for quantitative PCR.

*3 PrimeScript II RTase shows powerful extension ability even if the template RNA has a strong secondary structure. Therefore, it is generally recommended to perform the RT reaction at 42°C with this enzyme. However, in case of RT-PCR, if the reverse primer for PCR is used as a RT primer, nonspecific products may be amplified due to mispriming. In such a case, it is recommended to perform the RT reaction at 45 - 50°C.

The obtained reactant can be used for 2nd-strand synthesis reaction or as a template for PCR.

PrimeScript is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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