pRI 909 DNA

Code No. 3260

Size: 10 μ g

Conc.: 0.5 μ g/ μ l

* 2 years from date of receipt under proper storage conditions.

Description :

pRI 909 DNA is a binary vector which has a region of T-DNA for plant transformation. It is originated from Ri plasmid of *Rhizobium (Agrobacterium) rhizogenes*, lacking their vir region. pRI 909 DNA is also a shuttle vector, and replicates autonomously in *E. coli* and *Rhizobium (Agrobacterium)*. In *E. coli*, this vector is a high copy number plasmid because it has a replication origin same as that of pUC type plasmid (ColE1 ori), and it is maintained stably in *Rhizobium (Agrobacterium)* also with mutant type replication origin of Ri plasmid (Ri-ori). This vector has a kanamycin resistant gene, NPT III, as the selection marker for *E. coli* and *Rhizobium (Agrobacterium)*, and a mutant type kanamycin resistant gene, NPT II, as the selection marker for plant. This vector also has the same multi-cloning sites as pUC type plasmids, for insertion of the expression cassette of [promoter + target gene + terminator].

This is a vector for plant transformation using *Rhizobium (Agrobacterium)* which is available for binary vector method. It is possible to integrate target gene in plant chromosome stably using this vector because its cloning site is located at the position closer to Right Border (RB) of T-DNA than the selection marker (NPT II) for plant, so the target gene is not deleted.

Form : 10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage : -20°C


Preparation : Purification by column treatment

Base pairs : 9,168 bp



Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Cloning site :

cloning site


Usage :

This is a vector for transformation of plants such as Arabidopsis, Tomato, Tobacco, Rice using *Rhizobium (Agrobacterium)* which is available for binary vector method.

References :

- 1) R Nishiguchi, M Takanami, and A Oka. *Molecular and General Genetics*. (1987) **206**: 1-8.
- 2) G Ooms, J J Hooykaas, R J M V Veen, P V Beelen, T J G Regensburg-Tuink, and R A Schilperoort. *Plasmid*. (1982) **7**: 15-29.
- 3) T Ohba, Y Yoshioka, C Machida, Y Machida. *The Plant Journal*. (1995) **7**: 157-164.

Vector map of pRI 909 DNA :

ColE1 ori: Replication origin in *E. coli*

Ri ori: Replication origin in *Rhizobium (Agrobacterium)*

RB, LB: T-DNA border sequences

NOS-promoter, NOS-terminator :

Promoter or Terminator for expression of selectable marker gene in plant

NPT II: Selectable marker gene in plant

NPT III: Selectable marker gene in *E. coli* or *Rhizobium (Agrobacterium)*

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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