# T-Vector pMD19 (Simple)

**Code No. 3271**      **Size:**      **1  $\mu$ g**  
**Conc.:**      **50 ng/ $\mu$ l**

\* 2 years from date of receipt under proper storage conditions.

## Description :

T-Vector pMD19 (Simple) is a linearized vector with a single 3'-terminal thymidine at both ends. The T-overhang ends at the cloning site improve the efficiency of ligation of PCR products which contain A-overhangs at 3'-ends.


This vector have been deleted the multiple cloning sites in the *lacZ* gene sequence, but the  $\beta$ -galactosidase activity is not disrupted. Therefore, the transformant colony having a recombinant plasmid can also be identified by blue/white screening. Due to lack of multiple cloning sites on the vector, the inserted PCR product cannot be recut using restriction enzymes. When the insert DNA fragment is recut from the vector, the PCR primers should be added restriction site to these 5'-end for restriction enzyme digestion and subcloning.

T-Vector pMD19 (Simple) is supplied with a Control Insert (500 bp) DNA for positive control reaction.

## Components :

T-Vector pMD19 (Simple) (50 ng/ $\mu$ l)      20  $\mu$ l  
 Control Insert (50 ng/ $\mu$ l)      10  $\mu$ l

**Storage:**      -20°C

This product may be covered by licenses. Visit the product page for this product on our website to get the most updated information.

## Experimental Example :

1. Combine the following components in a microcentrifuge tube (total volume 5  $\mu$ l).

T-Vector pMD19 (Simple)	1 $\mu$ l
PCR product or Control Insert	1 $\mu$ l (0.1 - 0.3 pmol)*
Sterile purified water	3 $\mu$ l

2. Add one volume (5  $\mu$ l) of DNA Ligation Kit <Mighty Mix> (Cat. #6023) to the DNA solution and mix thoroughly.

3. Incubate the reactions for 30 minutes at 16°C.

**Notes :** (1) The ligation reaction can be finished in only 5 minutes, but the ligation efficiency may be slightly lower.

(2) Please extend the ligation reaction time to a few more hours if the fragments are over 2 kb.

4. Add the reaction solution (10  $\mu$ l) to 100  $\mu$ l of *E. coli* competent cells. Gently flick the microcentrifuge tube to mix and then place it on ice bath for 30 minutes.

5. Heat-shock the cells for 45 seconds in a water bath at 42°C. Immediately return the tube to ice bath for 1 minute.

6. Add 890  $\mu$ l of SOC Medium to the tube and incubate for 1 hour at 37°C with shaking.

7. Plate 100  $\mu$ l of transformation culture onto LB/ampicillin/IPTG/X-Gal plates, then incubate overnight at 37°C.

8. Select the white colonies, and confirm by PCR amplification.

\* The Control Insert DNA is a 500 bp PCR fragment with a single 3'-A overhang at both ends. 1  $\mu$ l of Control Insert DNA (50 ng) is approximately equal to 0.15 pmol. The ligation reaction has been optimized using a 1 : 2 to 10 molar ratio of the vectors to insert DNA.

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

```

lacZ →
GAGCGGA TAACAATTC ACACAGGAAA CAGCTATGAC CATGATTACG CCAAGTTTGC
CTCGCCT ATTGTTAAAG TGTGTCTTT GTGCATACTG GACTAATGC GGTCAAACG
    
```

```

EcoR V
ACGCCTGCCG TTCGAC GATT3' (cloned insert) ATC TCTGGAAG ATCCGCGCGT
TGGCGACGGC AAGCTG CTA 3'TTAG AGACCTTC TAGGCGCGCA
    
```

```

ACCGAGTTCT AATTCACTGG CCGTCGTTTT ACAACGTCGT GACTGGGAAA ACCCTGGCGT
TGCTCAAGA TTAAGTGACC GGCAGCAAAA TGTTGCAGCA CTGACCTTT TGGGACCGCA
    
```

M13 primer M4 ←

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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