# pUC118 *Pst* I/BAP

**Code No. 3323**

**Size:** 5  $\mu$ g

**Conc.:** 0.1  $\mu$ g/ $\mu$ l

\* 2 years from date of receipt under proper storage conditions.

## Description :

This DNA is a useful cloning vector, prepared by cleavage with *Pst* I which is a unique site in multi-cloning sites of pUC118, then dephosphorylated with alkaline phosphatase purified from *E. coli* (BAP). This DNA is able to prevent self-ligation of vectors, and also able to decrease blank value of transforming cells. Thus, this DNA is available to cloning without other treatment.



**Form :** 10 mM Tris-HCl, pH 8.0  
1 mM EDTA

**Storage:** -20°C

## Preparation :

Plasmid DNA pUC118 prepared by CsCl-EtBr ultracentrifugation was digested with *Pst* I, followed by dephosphorylation with Alkaline Phosphatase (BAP).

## Cloning site of pUC118 :

**Base pairs :** 3,162 bp  
(Calculated from the construction methods used by Messing, *et al.*)

**GenBank :** Entry Name : CVU07649  
Accession No. : U07649

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Application example :

pUC118 DNA <i>Pst</i> I/BAP	50 ng (25 fmol)
DNA fragment/ <i>Pst</i> I	25 - 250 fmol
TE buffer	up to 5 $\mu$ l
— DNA Ligation Kit <Mighty Mix> (Cat. #6023) Ligation Mix 5 $\mu$ l	
— 16°C, 30 min	
Transformation	

## Reference :

Vieira J and Messing J. *Methods in Enzymology*. (1987) **153**: 3-11.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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