α 1,2-L-Fucosidase

Code No. 4451

Size:

2 U

[Note]

The Specification of this product has been changed.

(Lot.N 104AA -)

- Form is changed to solution from lyophilized powder.
- The storage temperature is changed to -20°C from 4°C.

Systematic name: α -L-Fucoside fucosyltransferase

Enzyme code: 3.2.1.51

Source: *Corynebacterium* sp.

Intended usage: Specific cleavage of α 1,2-L-fucoside bonds

Form:

Solution in 50 mM potassium phosphate pH 8.0, 0.5 mM EDTA and 0.025% sodium azide.

Storage: -20°C

Avoid repeated Freeze-thawed. If necessary, aliquot into separate tubes and store at -20°C.

Concentration: 20 U/ml

Substrate specificity:

- This enzyme releases L-fucose from 2'- α -L-fucosyl lactose (Fuc α 1-2Gal β 1-4Glc) and from porcine gastric mucin oligosaccharides.
- This enzyme also releases L-fucose from a biantennary *N*-acetylglucosamin type sugar chain with α 1,6-fucoside bond (PA-Sugar Chain 009; Cat. #4109), but at a much lower reaction rate (Table 1).

Definition of activity:

One unit is the amount of enzyme required to hydrolyze 1 μ mol of *p*-nitrophenyl- α -L-fucoside (*p*-NP-Fucose) within 1 minute at 37°C at pH 8.5.

Assay of activity:

Incubated 0.7 mM *p*-nitrophenyl- α -L-fucoside as a substrate in 90 mM PBS, pH 8.5, at 37°C for 10 min. Stop the reaction by adding 0.25 M sodium carbonate. Calculate the enzyme activity on the basis of the absorbance at 405 nm of the reactant.

Properties:

Molecular weight	: 43,000
Michaelis constant	: $K_m=4.7 \times 10^{-5}$ M (<i>p</i> -NP- α -L-fucopyranoside)
Inhibitors	: Ag^+ , Hg^{2+} and Cu^{2+}
Optimum pH	: pH 8.5 (37°C 10 min)
Optimum temperature	: 34°C (pH 8.5, 10 min)
pH stability	: pH 5.5 - 9.5 (25°C 60 min)
Thermal stability	: $\leq 26^\circ C$ (pH 8.0, 60 min)

Purity:

- Protease Contamination:
No protease activity is detected after incubating 2 mU of α 1,2-L-Fucosidase with 4 nmol of oxidized insulin B chain for 16 hours at 37°C in 0.15 ml of 100 mM Tris-HCl buffer at pH 8.0.
- Endoglycosidase Contamination:
No endoglycosidase activity is detected after incubating 50 mU of α 1,2-L-Fucosidase with 100 pmol of pyridylamino oligomannose type sugar chain (M6-PA) for 16 hours at 37°C in 10 μ l of 100 mM sodium acetate buffer at pH 5.0.
- Exoglycosidase Contamination:
No detectable α -galactosidase, β -galactosidase, α -mannosidase, β -*N*-acetylglucosaminidase and β -xylosidase activities when 20 mU of α 1,2-L-Fucosidase is incubated with 500 nmol of *p*-nitrophenyl glycoside for 16 hours at 37°C in 0.5 ml of 100 mM sodium acetate buffer at pH 5.0.
- α 1,3/4-L-Fucosidase Contamination:
No α 1,3/4-L-Fucosidase activity is detected after incubating 600 mU of α 1,2-L-Fucosidase with 20 pmol of pyridylamino lacto-*M*-fucopentaosell or III for 63 hours at 37°C in 10 μ l of 30 mM potassium phosphate buffer at pH 7.0. (Table 1)

References:

- 1) Sano M, Sakai T, Hayakawa K, and Kato I. *Glycoconjugate J.* (1991) **8**: 271.
- 2) Yamashita K, Totani K, Iwaki Y, Kuroki M, Matsuoka Y, Endo T, and Kobata A. *J Biol Chem.* (1989) **264**: 17873.
- 3) Amano J, Straehl P, Berger E, Kochibe N, and Kobata A. *J Biol Chem.* (1991) **266**: 11461.
- 4) Matsui T, Titani K, and Mizuochi T. *J Biol Chem.* (1992) **267**: 8723.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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