SCDase

(Sphingolipid ceramide *N*-deacylase)

Code No. 4462

Size: 250 mU/50 μ l

Description:

SCDase (Sphingolipid ceramide *N*-deacylase) from *Pseudomonas* sp. hydrolyzes the *N*-acyl linkage between fatty acids and sphingosine bases in ceramide of various sphingolipid.¹⁾ The enzyme catalyzes the reverse reaction and transacylation, too.²⁾

The enzyme well acts on various acidic and neutral glycosphingolipids and sphingomyelin. But it does not well act on ceramide.

Source: *Pseudomonas* sp.

Form:

Solution in 50 mM sodium acetate pH6.0 containing 0.1% Lubrol PX

Storage:

Store at or below -20°C until use.

Once it thawed, it is recommended to dispense the solution and store at -20°C. Repeated freeze-thaw cycles should be avoided.

Properties:

Optimum pH: Hydrolysis pH5.0 - 6.0 (0.8% Triton X-100)

Reverse reaction pH7.0 (0.1% Triton X-100)

Inhibitors: Hg²⁺, Zn²⁺ and Cu²⁺.

The reaction is not inhibited by Mn²⁺, Ca²⁺, Mg²⁺ and EDTA.

Unit Definition:

One unit (U) of enzyme was defined as the amount needed to catalyze the hydrolysis of 1 μ mol of GM1 per minute.

Assay Method:

After 10 μ l of the reaction mixture including this enzyme, 2 nmol of GM1, 0.8% Triton X-100 and 50 mM acetate buffer, pH6.0 is incubated at 37°C for 15 min., it is analyzed by thin layer chromatography. GM1 and lyso-GM1 are visualized by spraying the TLC plates with orsinol-H₂SO₄ reagent and scanning them with image analyzer.

Purity:

It was confirmed that this enzyme is free from the following exoglycosidase activities:

α -galactosidase

β -galactosidase

α -*N*-acetylgalactosaminidase

β -*N*-acetylgalactosaminidase

β -*N*-acetylglucosaminidase

α -mannosidase

α -fucosidase

sialidase

The preparation did not show protease and sphingomyelinase activities.

Assay Method

1) Exoglycosidase contamination test:

No exoglycosidase activity was detected by HPLC analysis when 20 mU SCDase is incubated with 20 pmol of a pyridylamino oligosaccharide for 16 hrs. at 37°C, in 10 μ l of acetate buffer, pH6.0.

2) Protease contamination test:

No protease activity was detected after incubation of 50 mU SCDase with 40 μ g of casein labeled with resorufin for 16 hrs. at 37°C, in 40 μ l of acetate buffer, pH6.0.

3) Sphingomyelinase contamination test:

No sphingomyelinase activity was detected after incubation of 20 mU SCDase with 15 nmol of TNPAL-sphingomyelin for 16 hrs. at 37°C, in 55 μ l of acetate buffer, pH6.0.

References:

- 1) Ito M, Kurita T, and Kita K. *J Biol Chem.* (1995) **270**: 24370-24374.
- 2) Sueyoshi N, Izu H, and Ito M. *J Lipid Res.* (1997) **38**: 1923-1927.
- 3) Mitsutake S, Kita K, Okino N, and Ito M. *Anal Biochem.* (1997) **247**: 52-57.
- 4) Tani M, Kita K, Komori H, Nakagawa T, and Ito M. *Anal Biochem.* (1998) **263**: 183-188.
- 5) Mitsutake S, Kita K, Nakagawa T, and Ito M. *J Biochem.* (1998) **123**: 859-863.
- 6) Izu H, Kurita K, Sano M, Kato I, and Ito M. XIX International Carbohydrate Symposium Abstract Book. (1998) Aug: DP066.
- 7) Kurita K, Izu H, Sano M, Ito M, and Kato I. *J Lipid Res.* (2000) **41**: 846-851.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

SCDase

(Sphingolipid ceramide *N*-deacylase)

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