Template Vector (*BspQ I*) for T7 mRNA Synthesis

Code No. 6146

Size:

10 μ l

Description:

Template Vector (*BspQ I*) for T7 mRNA Synthesis is used for the construction of template plasmid for *in vitro* transcription (IVT) using cap analogs [CleanCap Reagent AG or CleanCap Reagent AG (3' OMe)] from TriLink BioTechnologies. The pre-linearized vector contains a T7 promoter, transcription start sequence (AGG), 5'-UTR (untranslated region), 3'-UTR, and a 141-base Poly(A) sequence, simply allowing the construction of a IVT template plasmid just by cloning a desired gene (coding sequence; CDS) with In-Fusion® (Figure. 1, 2). Followed by a digestion with Type IIS restriction enzyme, *BspQ I*, the resulting linearized template plasmid makes it possible to provide a scarless mRNA after the Poly(A) sequence, to be used for human, murine, or other mammalian cells. The linearized template plasmid can then be used to produce the desired mRNA by performing high-yield *in vitro* transcription using Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144).

<In-Fusion cloning>


Fig. 1. IVT template plasmid construction by In-Fusion cloning


Fig. 2. Transcription start sequence (AGG) of IVT template when CleanCap Reagent AG is used

Components:

Linearized Template Vector (<i>BspQ I</i>) (50 ng/ μ l)	10 μ l
FLuc Control Fragment (<i>BspQ I</i>) (100 ng/ μ l)	10 μ l

Storage: -20°C

DNA Length:

2,962 bp [Linearized Template Vector (*BspQ I*)]
1,683 bp [FLuc Control Fragment (*BspQ I*)]

* For more information on the sequences, please see our website.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Usage:

Preparation of linearized template plasmid for mRNA with scarless sequence after Poly(A)

Precautions for Use:

1. A cloning reagent for In-Fusion cloning is required for the IVT template plasmid construction. Please purchase In-Fusion Snap Assembly Master Mix separately (Cat. #638943/638944/638947 - 638949).
2. Please refer to the user manual of Cloning Kit for mRNA Template (Cat. #6143) for more detailed information on the design of template plasmid and procedures. Confirm that the *BspQ I* restriction site used for linearizing the IVT template plasmid is not present in the CDS of the desired gene. If the restriction site *BspQ I* is present in the CDS, change the DNA sequence by changing the codon in the recognition sequence while keeping the amino acid sequence the same (e.g., switch from alanine codon "GCU" to "GCC").
3. RNase contamination of the double-stranded DNA template, reagents, tubes, micropipette tips, or other materials used in the reaction can significantly decrease or digest RNA obtained with the kit. Use dedicated tubes and micropipette tips in the reaction and wear new disposable gloves to prevent RNase contamination.

References:

- 1) Karikó, K. *et al.* Incorporation of pseudouridine into mRNA yields superior nonimmunogenic vector with increased translational capacity and biological stability. *Mol Ther.* (2008) **16**: 1833-1840.
- 2) Vaidyanathan, S. *et al.* Uridine Depletion and Chemical Modification Increase Cas9 mRNA Activity and Reduce Immunogenicity without HPLC Purification. *Mol Ther Nucleic Acids.* (2018) **12**: 530-542.
- 3) Xia, X. Detailed Dissection and Critical Evaluation of the Pfizer/BioNTech and Moderna mRNA Vaccines. *Vaccines (Basel).* (2021) **9**: 734.
- 4) Schenborn, E. T. and Mierendorf, R. C. A novel transcription property of SP6 and T7 RNA polymerases: dependence on template structure. *Nucleic Acids Res.* (1985) **13**: 6223-6236.

Related Products:

In-Fusion® Snap Assembly Master Mix

(Cat. #638943/638944/638947 - 638949)

BspQ I (Cat. #1227A/B)

PrimeSTAR® Max DNA Polymerase (Cat. #R045A/B)

TaKaRa Ex Premier™ DNA Polymerase (Cat. #RR370S/A/B)

TaKaRa Ex Premier™ DNA Polymerase Dye plus (Cat. #RR371S/A/B)

NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/50/.250)*

NucleoSpin Plasmid (Cat. #740588.10/50/.250)*

Takara IVTpro™ mRNA Synthesis System (Cat. #6141)

Cloning Kit for mRNA Template (Cat. #6143)

Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)

RNase-free Water (Cat. #9012)

T7 RNA Polymerase ver.2.0 (Cat. #2541A)

T7 RNA Polymerase, HQ (Cat. #2542A)

Pyrophosphatase (inorganic) (Cat. #2450A/B)

*: Not available in all geographic locations. Check for availability in your area.

PrimeSTAR is a registered trademark of Takara Bio Inc.

In-Fusion is a registered trademark of Takara Bio USA, Inc.

IVTpro and Ex Premier are trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Template Vector (*BspQ* I) for T7 mRNA Synthesis

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− 20°C

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2,962 bp [Linearized Template Vector (*BspQ* I)]
1,683 bp [FLuc Control Fragment (*BspQ* I)]

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