Protein disulfide-isomerase (PDI)


Code No. 7318

Size:

1 mg

Description:

This enzyme accelerates the exchange reactions between disulfide bonds in proteins in the presence of appropriate oxidizing or reducing reagents.

Source: Bovine liver

Form: Lyophilized white powder
(containing the amount equivalent to 200 μ l of 50 mM sodium phosphate buffer, pH 7.5)
Dissolve the content in 200 μ l of sterile purified water.

Storage: -20°C

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Definition of Activity:

One unit of PDI is defined as the amount required to recover one RNase A unit from reduced bovine RNase A, at 25°C, pH 7.5 for 15 minutes. One RNase A unit is defined as the amount required for hydrolysis of cCMP (cytidine 2':3'-cyclic monophosphate) that results in an increase of 1×10^{-3} absorbance unit at 284 nm per minute.

Activity: > 300 U/mg protein

Properties:

Molecular weight: 107,000 (homodimer)
Isoelectric point: 4.2
Optimum pH: pH 7.0 - 9.0

Assay of activity:

<Principle>

Reduced bovine RNase A is treated with PDI in the presence of cCMP (cytidine 2':3'-cyclic monophosphate). Reactivation of RNase A is monitored by the absorbance at 284 nm.

<Protocol>

Reagents

- | | | |
|----|-----------|--|
| A. | 1 mg/ml | Reduced RNase A |
| B. | 0.2 M | Sodium phosphate buffer, pH 7.5 |
| C. | 20 mM | GSH (6.15 mg/ml) |
| D. | 2 mM | GSSG (1.23 mg/ml) |
| E. | 0.1 mg/ml | cytidine 2':3'-cyclic monophosphate/0.1 M MOPS, pH 7.0 |

Procedure

1. Mix 40 μ l of A, 100 μ l of B, 20 μ l of C and 20 μ l of D, then add 20 μ l of PDI (< 1 mg/ml) to the mixture.
2. Incubate at 25°C for 15 min.
3. Add 2.8 ml of E.
4. Measure the absorbance at 284 nm.
5. Incubate at 25°C for 3 min.
6. Measure the absorbance at 284 nm.

Calculations

$$\text{Activity (U/ml)} = \frac{\Delta A_{284} \times 10^3 \times \text{DF}}{\text{Rt} \times \text{Ev}}$$

ΔA_{284} : $\Delta A_{284} (\text{complete}) - \Delta A_{284} (-\text{PDI}) - \Delta A_{284} (-\text{reduced RNase A})$
DF: dilution factor
Rt: RNase A reaction time (3 min)
Ev: Amount of enzyme (0.02 ml)

Application:

Useful for refolding of proteins by exchange reactions between disulfide bonds.

References:

1. Goldberger R F, Epstein C J, and Anfinsen C B. *J Biol Chem.* (1964) **239**: 1406.
2. Lambert N and Freedman R B. *Biochem J.* (1983) **213**: 225.
3. Hillson D A, Lambert N, and Freedman R B. *Methods in Enzymology.* (1984) **107**: 281.
4. Tang J G, Wang C C, and Tsou C L. *Biochem J.* (1988) **255**: 451.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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