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Improved Accurate Quantitative Analysis of Microbiome Using DNA Standard for 16S rRNA NGS analysis

Honami Miyakura, Yasunobu Terabayashi, Yuto Kitajima, Tomoya Uchiyama, Yoshitaka Kimura
TAKARA BIO INC.

Summary

High-throughput 16S rRNA gene amplicon analysis (16S rRNA microbiota analysis) is widely used for microbiota profiling. However, previous studies have been limited to relative quantification due to the lack of reliable internal standards, which are indispensable for absolute quantification of each species in the microbiota.


In this study, to construct widely applicable NGS-based absolute quantitative 16S rRNA microbiota analysis protocol, availability of the synthetic 16S rDNA set as internal standard was evaluated, in conjunction with a data analysis workflow. This set consists of 12 known artificial DNA sequences in different known proportions, and each standard possesses a unique artificial sequence not found in nature in the variable region (V1-V9) of 16S rDNA^[1]. To employ this set as a spike-in control, five required characteristics were evaluated: reliability, robustness with various regions, consistency, reproducibility, and versatility.

As a result, a widely applicable absolute quantitative NGS-based microbiota analysis protocol was established with using synthetic 16S rDNA standard set as spike-in control.

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

Conventional 16S rRNA microbiota analysis: relative abundance of each bacterial species.

Absolute quantification of 16S rRNA microbiota analysis with synthetic 16S rDNA standards enables:

- ✓ To observe the copy numbers of each bacterial 16S rRNA.
- ✓ To compare bacterial abundance among different experiment lots.
- ✓ To define detection limits based on copy numbers.
- ✓ To control the quality of analysis process.

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

*: ZymoBIOMICS : Microbial Community DNA Standard

Conclusion

Five characteristics of "Synthetic 16S rDNA standards set" required for spike-in control were evaluated.

1. Correlation of quantitative results with other techniques (ddPCR)
2. Reproducibility of quantification results from various regions (V1-V2, V3-V4, and V4)
3. Consistency of results of 16S rRNA microbiota analysis with and without standards
4. Reproducibility of results in repeated experiments
5. Versatility of spike-in for various samples

The result indicated broad availability of synthetic 16S rDNA standard set as spike-in controls for NGS-based 16S rRNA absolute quantitative microbiota analysis protocol. Absolute quantification with this protocol will lead to additional interpretations and insights to microbiota profiling.

Reference

1. Tourlousse D, et al. Synthetic spike-in standards for high-throughput 16S rRNA gene amplicon sequencing. *Nucleic Acids Res* . 2017 Feb 28; 45(4): e23.
2. Klindworth A, et al. Evaluation of general 16S ribosomal RNA gene PCR primers for classical and next-generation sequencing-based diversity studies. *Nucleic Acids Res* .2013 Jan 7, 41(1): e1.
3. Caporaso JG, et al. Ultra-high-throughput microbial community analysis on the Illumina HiSeq and MiSeq platforms. *ISME J* . 2012 Aug, 6(8): 1621-1624.