Sample Protector for RNA/DNA

Code No. 9750

Size: 100 ml

Description:

Sample Protector for RNA/DNA is a non-toxic stabilization solution for storage of mammalian or plant tissues, cultured cells, yeast, or bacteria prior to DNA or RNA extraction. It rapidly inactivates RNase and DNase, allowing these samples to remain stable at room temperature or 4°C for specific periods of time (Please refer to Protocol 2. Sample storage). When RNA or DNA extraction from fresh tissue or cell samples cannot be performed immediately after harvesting, the samples are often stored at -20°C, which can result in some RNA and DNA degradation. This product can effectively solve common problems with storage and shipment of fresh tissue and cell samples to be used for DNA or RNA preparation. In addition, if biological samples collected at various time points are immediately submerged in Sample Protector for RNA/DNA, this reagent can stabilize RNA and preserve gene expression profiles, enabling gene expression analysis. The effect of preservation is equal to that in -80°C refrigerator for specific period of time.

Storage: 4°C

This product will remain stable until the expiration date, but performance may diminish with longer-term storage.

Protocol:

1. Sample processing

(1) Mammalian and plant tissues

Cut the mammalian or plant tissue into 0.5-cm square pieces and add 0.5 - 1 ml of Sample Protector for RNA/DNA, mix upside down and immerse tissue under the surface as much as possible. Then store under appropriate conditions as described in Step 2. Sample storage.

(2) Cultured cells and leukocytes

Harvest the cells ($10^6 - 5 \times 10^6$, centrifuge <5,000 rpm, avoid cells damage as much as possible) according to the standard protocol, wash with 0.5 - 1 ml of PBS, centrifuge <5,000 rpm, (avoid cells damage as much as possible), discard supernatant, add 0.5 - 1 ml of Sample Protector for RNA/DNA, and completely suspend cells. Store under appropriate conditions as described in Step 2. Sample storage.

(3) Adherent cells

Treat cells with trypsin, harvest the cells ($10^6 - 5 \times 10^6$) by centrifuging <5,000 rpm, (avoid cells damage as much as possible) and discard the supernatant. Wash with 0.5 - 1 ml of PBS, centrifuge <5,000 rpm, (avoid cells damage as much as possible) and discard the supernatant. Add 0.5 - 1 ml of Sample Protector for RNA/DNA and completely suspend cells. Store under appropriate conditions as described in Step 2. Sample storage.

(4) Yeast

Harvest approximately 10^8 cells by centrifuging at 12,000 rpm for 2 min and discard the supernatant. Add 0.5 - 1 ml of Sample Protector for RNA/DNA and completely suspend the cells. Store under appropriate conditions as described in Step 2. Sample storage.

(5) Bacteria

Harvest approximately 10^9 cells by centrifuging at 12,000 rpm for 2 min and discarding the supernatant. Add 0.5 - 1 ml of Sample Protector for RNA/DNA, completely suspend the cells. Store under appropriate conditions as described in Step 2. Sample storage.

2. Sample storage

Sample can be stored in Sample Protector for RNA/DNA for up to 1 day at 37°C, 1 week at 25°C, 1 month at 4°C, or long-time storage* at -20°C.

* For long-time storage, incubate the samples in Sample Protector for RNA/DNA for 1 hour at room temperature then store at -20°C.

3. RNA or DNA extraction

- (1) Remove the cell or tissue sample for RNA or DNA extraction from the Sample Protector for RNA/DNA solution as follows: Remove tissue samples from the Sample Protector solution using sterile forceps, etc. Separate cultured cells from Sample Protector solution by centrifuging at >5,000 rpm for 5 min or 12,000 rpm for 1 min, as Sample Protector for RNA/DNA is the higher density.
- (2) Wash the tissues or cells once or twice with PBS or RNase-free H₂O, by centrifuging at appropriate centrifugal force.
- (3) Extract RNA or DNA from the sample using a conventional nucleic acid extraction kit.

Note:

1. Small tissue samples such as mouse kidney or spleen can be stored in Sample Protector for RNA/DNA directly without cutting into smaller pieces.
2. Sample Protector for RNA/DNA is a non-toxic reagent.
3. For long-time storage, incubate the samples in the Sample Protector for RNA/DNA for 1 hour at room temperature, then store at -20°C.
4. When the sample is stored in Sample Protector for RNA/DNA at -20°C, it may be in the state of solution or freezing.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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
Sample Protector for RNA/DNA


Experimental Examples:

1. Extraction of total RNA from HL-60 cells that were stored in Sample Protector for RNA/DNA.

- 1: Stored at 37°C for 1 day after adding Sample Protector
2: Stored at 25°C for 1 week after adding Sample Protector
3: Stored at 37°C for 1 day without Sample Protector

2. Extraction of genomic DNA from rat spleen that was stored in Sample Protector.

- 1: Stored at 37°C for 1 day after adding Sample Protector
2: Stored at 25°C for 1 week after adding Sample Protector
3: Stored at 37°C for 1 day without Sample Protector

Q & A:

- Q1: Which samples are suitable for storage in the Sample Protector for RNA/DNA?
A1: Sample Protector for RNA/DNA is suitable for storage of fresh samples. For frozen samples, add Sample Protector for RNA/DNA before melting of the samples (that is, just be taken out from the freezer or liquid nitrogen), vortex to make samples thawing, then store at $\geq 4^{\circ}\text{C}$; If further cryopreservation is required, the sample will be placed in freezing conditions immediately after adding Sample Protector for RNA/DNA.
Note: If adding this product after samples melt, the storage effect will become worse.
- Q2: Is Sample Protector for RNA/DNA suitable for storage of whole blood?
A2: No. Because a lot of insoluble matter will appear after adding this product in the whole blood, it can't store the whole blood. Lymphocyte separation solution can be used to separate leukocytes, and then store with this product.
- Q3: Can samples be stored in the Sample Protector for RNA/DNA at -80°C ?
A3: Under the condition of -80°C , the samples can be frozen directly. It is not recommended to add Sample Protector for RNA/DNA.

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