# Tks Gflex™ DNA Polymerase

**Code No. R060A**

**Size: 250 units  
(200 rxn)**

**Shipping at – 20°C  
Store at – 20°C**

## Components :

<b>Tks Gflex DNA Polymerase</b>	<b>250 U</b>
<b>2X Gflex PCR Buffer (Mg<sup>2+</sup>, dNTP plus)*</b>	<b>1 ml x 5</b>
<b>* : Mg<sup>2+</sup> conc. (2X)</b>	<b>2 mM</b>
<b>dNTP conc. (2X)</b>	<b>ea. 400 μM</b>

## Lot No.

**Conc. : 1.25 units/μl**

**Volume : 200 μl**

## Expiration Date :

### Description :

Tks Gflex DNA Polymerase is a PCR enzyme improved based on a high fidelity DNA polymerase derived from *Thermococcus* sp. This enzyme possesses high fidelity of a family B DNA polymerase and an excellent extension activity by suppressing its non-specific binding to template DNA. The Takara Bio's original strong elongation factor and the newly developed accelerator for specific priming (patent pending) contribute to high-speed amplification with extremely high specificity from templates in wide range of nucleic acid concentration. Therefore this polymerase makes PCR amplification succeeded in most cases with such targets much longer and that contain GC- or AT-rich sequences difficult to amplify by general DNA polymerases.

The elements which absorb the substances inhibiting PCR or accelerate amplification in the reaction buffer enable PCR using crude samples to be in high efficiency.

This polymerase is for Hot Start PCR utilizing a monoclonal antibody which suppresses DNA polymerase activity and 3' -5' exonuclease activity at room temperature.

<b>Storage Buffer:</b>	50 mM Tris-HCl, pH8.2 at 4°C
	100 mM NaCl
	0.1 mM EDTA
	1 mM DTT
	0.1% Tween 20
	0.1% Nonidet P-40
	50% Glycerol

### Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

### Purity :

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 1 μg of supercoiled pBR322 DNA, 1 μg of λ DNA or 1 μg of λ-Hind III digest with 15 units of this enzyme for 1 hour at 74°C.

### Application:

DNA amplification by Polymerase Chain Reaction (PCR).

This enzyme is recommended for :

- PCR amplification with GC-, or AT-rich templates.
- PCR amplification of long targets as more than 10 kb.
- PCR amplification with crude extracts.
- PCR amplification with cDNA containing a large amount of nucleic acids.

**PCR products :** Most of PCR products obtained using Tks Gflex DNA Polymerase will possess blunt-ends. The obtained PCR products can be directly cloned into blunt-end vectors. (If necessary, phosphorylate PCR products before cloning.)

**Electrophoresis :** TAE buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using Tks Gflex DNA Polymerase. Use of TBE buffer may result in DNA band patterns which are enlarged at the gel bottom.

**PCR test :** Good performance of this product was confirmed by PCR using human genomic DNA as a template (amplified fragments; 0.5, 4, 30 kb).

### General reaction mixture for PCR (total 50 μl) :

Tks Gflex DNA Polymerase (1.25 units/μl)	1 μl
2X Gflex PCR Buffer (Mg <sup>2+</sup> , dNTP plus)	25 μl
Template	< 500 ng
Primer 1	0.2 - 0.3 μM (final conc.)
Primer 2	0.2 - 0.3 μM (final conc.)
Sterilized distilled water	up to 50 μl

### PCR conditions (an example):

98°C	10 sec.	} 30 cycles	or	98°C	10 sec.	} 30 cycles
60°C	15 sec.		68°C	30 sec./kb		
68°C	30 sec./kb					

**(Note)** Annealing temperature depends on the Tm of primers for PCR. The initial denaturation of templates at 94°C for 1 min. is recommended for long or GC-rich targets. The extension time of 1 min./kb is recommended for crude samples.

NOTICE TO PURCHASER: LIMITED LICENSE

### [P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

### [L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287, and corresponding patents in other countries.

### [M54] PrimeSTAR HS DNA Polymerase

This product is covered by the claims of U.S. Patent Nos. 7,704,713 and its foreign counterparts.

### Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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