# TaKaRa Ex Taq<sup>®</sup>

Code No. RR001A  
Size: 250 units

Shipping at  $-20^{\circ}\text{C}$   
Store at  $-20^{\circ}\text{C}$

## Supplied Reagents :

10X Ex Taq Buffer (20 mM Mg<sup>2+</sup> plus) 1 ml  
dNTP Mixture (2.5 mM each) 800  $\mu\text{l}$

Conc. : 5 units/ $\mu\text{l}$   
Volume : 50  $\mu\text{l}$

## Storage buffer :

20 mM Tris-HCl (pH8.0)  
100 mM KCl  
0.1 mM EDTA  
1 mM DTT  
0.5% Tween 20  
0.5% Nonidet P-40  
50% Glycerol

## Supplied dNTP Mixture :

dNTP Mixture is ready for use in PCR without dilution.

Concentration : 2.5 mM of each dNTP  
Form : Dissolved in water (sodium salts), pH7 - 9  
Purity :  $\geq 98\%$  for each dNTP

## Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

## Reaction mixture for unit definition :

25 mM TAPS (pH9.3 at 25°C)  
50 mM KCl  
2 mM MgCl<sub>2</sub>  
0.1 mM DTT  
200  $\mu\text{M}$  each dATP · dGTP · dCTP  
100  $\mu\text{M}$  [<sup>3</sup>H]-dTTP  
0.25 mg/ml activated salmon sperm DNA

## Purity :

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 0.6  $\mu\text{g}$  of supercoiled pBR322 DNA, 0.6  $\mu\text{g}$  of  $\lambda$  DNA or 0.6  $\mu\text{g}$  of  $\lambda$ -Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

**Applications :** For DNA amplification by polymerase chain reaction (PCR)

## PCR products :

Most PCR products amplified with TaKaRa Ex Taq have one A added at the 3'-terminus. Thus, the PCR product can be used directly for cloning into a T-vector. Additionally, it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

## PCR test :

Good performance of DNA amplification by PCR was confirmed using  $\lambda$  DNA as the template (amplified fragment: 20 kb).

Good performance of DNA amplification of a  $\beta$ -globin gene fragment by PCR was also confirmed using human genomic DNA as the template (amplified fragment: 17.5 kb).

## General reaction mixture for PCR (total 50 $\mu\text{l}$ ) :

TaKaRa Ex Taq (5 units/ $\mu\text{l}$ )	0.25 $\mu\text{l}$
10X Ex Taq Buffer	5 $\mu\text{l}$
dNTP Mixture (2.5 mM each)	4 $\mu\text{l}$
Template	< 500 ng
Primer 1	0.2 - 1.0 $\mu\text{M}$ (final conc.)
Primer 2	0.2 - 1.0 $\mu\text{M}$ (final conc.)
Sterilized distilled water	up to 50 $\mu\text{l}$

## Example of PCR conditions :

When amplifying a 1 kb DNA fragment

98°C 10 sec.	] 30 cycles	or	98°C 10 sec.	] 30 cycles
55°C 30 sec.		68°C 1 min.		
72°C 1 min.				

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. We recommend denaturing for 5 - 10 sec. at 98°C or 20 - 30 sec. at 94°C.

## < Cool Start Method >

The Cool Start Method provides more accurate amplification and minimizes amplification of nonspecific bands. This is a simple method that does not require specialized enzymes or additional reagents.

## Protocol for Cool Start Method

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice.\*1,2
  - \* 1 : The order in which reagents are added does not influence results.
  - \* 2 : Results will not be affected by leaving the mixture on ice for 30 min. before thermal cycling.
- 3) Program a thermal cycler with the desired cycling conditions.\*3
  - \* 3 : PCR conditions do not need to be altered for the Cool Start Method.
- 4) Set the tubes in a thermal cycler and start thermal cycling immediately.

TaKaRa Ex Taq is a registered trademark of TAKARA BIO INC.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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