

NucleoSpin® RNA – Support Protocol – RNA from PAXgene® tubes

Protocol Details

Application	RNA isolation from PAXgene® tubes	
Kit	NucleoSpin® RNA	
REF	740955 (.10/ .50/ .250)	
Protocol name	Support Protocol – RNA from PAXgene® tubes	Rev01



Protocol Steps

Steps	Procedure
1. Harvest sample	Centrifuge the PAXgene® Blood RNA Tube for 10 min at 3000 – 5000 x g using a swing-out rotor.
2. Wash pellet	Remove the supernatant by decanting or pipetting. Add 4 mL RNase-free water to the pellet, and close the tube using a fresh secondary Hemogard closure. Thoroughly resuspend the pellet by vortexing, and centrifuge for 10 min at 3000 – 5000 x g. Remove and discard the entire supernatant. Gently decant the supernatant, and blot on a paper towel. Thoroughly resuspend the pellet in 300 µL RNase-free water by pipetting. Do not vortex.
3. Lyse sample	Pipette the sample into a 1.5 ml microcentrifuge tube. Add 300 µL RA1 supplemented with 3 µL β-Mercaptoethanol, and 40 µL Proteinase K . Mix by vortexing, and incubate for 20 min at 55 °C using a shaker incubator, heating block, or water bath.
4. Filtrate sample	Clear the lysate by filtration through NucleoSpin® Filter (violet rings). Place the NucleoSpin Filter in a Collection Tube (2ml), apply the mixture and centrifuge for 1 min at 11,000 x g.
5. Adjust RNA binding conditions	Discard the NucleoSpin® Filter and transfer flow-through (cleared lysate) to a new reaction tube. Add 300 µL ethanol (96 %) to the homogenized lysate and mix by pipetting up and down (5 – 10 times).
6. Transfer lysates	Load the lysate to the NucleoSpin® RNA column (light blue ring). Proceed with step 5 of the standard protocol ('Bind RNA').

www.mn-net.com



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com

Reagents, consumables, and equipment to be supplied by user

Equipment	Content	REF
Reagents		
Proteinase K	2x50 mg Proteinase K 8 ml Proteinase Buffer PB	740506

Equipment

Water bath or
heating block

Trademarks

NucleoSpin® is a trademark of MACHEREY-NAGEL GmbH & Co. KG

PAXgene® is a trademark of PreAnalytiX

