

## Supplementary protocol

---

### NucleoSpin® Plasmid – isolation of M13 DNA (Rev. 01)

*This protocol is only a supplement to the kit's general user manual. Please refer to the kit manual for more detailed information regarding safety instructions, product-specific disclaimers, and especially preparations needed before starting the procedure. The latest version of the user manual is available at [www.mn-net.com/usermanuals](http://www.mn-net.com/usermanuals) or can be requested from our technical service ([tech-bio@mn-net.com](mailto:tech-bio@mn-net.com)). Material safety data sheets (MSDS) can be downloaded from [www.mn-net.com/MSDS](http://www.mn-net.com/MSDS).*

#### Additional reagent needed:

- Buffer NT2 (REF 740597)
  - Glacial acetic acid
- 

#### 1 Cultivate and harvest bacterial cells

Grow a M13-infected *E. coli* culture and centrifuge bacterial cells at **4,000 x g** for **10 min** at **4 °C**.

---

#### 2 Cell lysis

Transfer the supernatant to a new microcentrifuge tube. For each **1 mL supernatant** (containing suspended phage particles), add **10 µL glacial acetic acid**. Mix by inverting the tube 6–8 times. Incubate at room temperature for **2 min**.

Place a NucleoSpin® Plasmid Column into a Collection Tube (2 mL) and load sample. Centrifuge for **1 min** at **11,000 x g** and discard flow-through.

*Maximal loading volume of a NucleoSpin® Plasmid Column is 700 µL. If larger volumes are to be processed, load samples in successive steps. Do not load the column more than 3 times.*

Place the NucleoSpin® Plasmid Column back into the collection tube and add **700 µL** of **Buffer NT2**. Centrifuge for **1 min** at **11,000 x g** and discard flow-through.

Place the NucleoSpin® Plasmid Column back into the collection tube, add 700 µL of Buffer NT2, and incubate for 1 min at room temperature. Centrifuge for 1 min at 11,000 x g and discard flow-through.

---

#### 3 Isolate M13 DNA

Continue with step 5 ('Wash silica membrane') of the NucleoSpin® Plasmid standard protocol.

---