

## Accurate aneuploidy detection with the Embgenix PGT-A Kit


- The Embgenix PGT-A Kit is an end-to-end solution for next-generation sequencing (NGS)-based preimplantation genetic testing of embryos for aneuploidies (PGT-A)
- The kit integrates [PicoPLEX whole-genome amplification \(WGA\)](#), Illumina library preparation, and automated analysis software to streamline the assessment of trophoctoderm (TE) biopsy and genomic DNA (gDNA) samples for whole chromosome and segmental aneuploidies and mosaicisms
- The kit features a streamlined workflow with less than 2 hours of hands-on time and a single dilution factor following WGA; unique dual indexes are provided in a ready-to-use plate format to enable preparation of Illumina libraries and multiplexing of up to 96 samples
- Cloud-based Embgenix Analysis Software employs a proprietary algorithm for automated assessment of copy number variations (CNVs) from sequencing data generated with the Embgenix PGT-A Kit
- The Embgenix PGT-A workflow enables fast and accurate detection of segmental aneuploidies  $\geq 10$  Mb and low (30–50%) and high (50–70%) mosaicisms

### Introduction ^

The aim of this study was for multiple laboratories—Takara Bio's in-house laboratory and independent clinical laboratories A, B, and C—to test the performance of the Embgenix PGT-A Kit in detecting copy number variations (CNVs), such as whole chromosome, segmental, and mosaic aneuploidies in TE biopsy and gDNA samples.


NGS is the method of choice for PGT-A of in vitro fertilization (IVF) embryos (Vendrell et al. 2017). The overall goal of PGT-A is to improve rates of embryo implantation and live birth by enabling assessment and prioritization of embryos for transfer on the basis of copy number status in addition to other established criteria (Kemper et al. 2019). With integrated workflows from sample preparation to analysis, NGS-based PGT-A offers a simpler, more scalable solution for aneuploidy screening as compared to preceding methods.

Leveraging Takara Bio's technologies for WGA and library preparation and a customized bioinformatics pipeline, Embgenix PGT-A is an end-to-end solution for detection of whole chromosome, segmental, and mosaic aneuploidies in human embryos. The Embgenix PGT-A Kit has been optimized to analyze DNA obtained from embryo biopsies or genomic DNA. The workflow starts with cell lysis and PicoPLEX WGA from template DNA followed by library preparation and incorporation of unique dual indexes (UDI). Following purification, quantification, and pooling, libraries are ready for multiplex sequencing on Illumina NGS platforms (**Figure 1**).




**Figure 1.** Embgenix PGT-A workflow. Starting from cells or gDNA, up to 96 Illumina-ready sequencing libraries can be generated in parallel in six to eight hours.

PicoPLEX WGA, used to generate enough DNA for NGS library preparation, is shown in **Figure 2**. PicoPLEX is the gold standard for PGT-A applications by virtue of its superior reproducibility, fidelity, and uniformity as compared to other WGA methods (Zhang et al. 2017). WGA products are enzymatically fragmented and then stem-loop adapters and UDI are incorporated to yield high-quality Illumina-ready libraries (**Figure 3**). Sequencing and analysis of the libraries enables accurate screening of all 24 chromosomes for whole- and sub-chromosomal CNVs.



**Figure 2.** Takara Bio's PicoPLEX WGA chemistry. WGA products are generated from gDNA extracted from TE biopsy samples via cell lysis. Pre-amplification and amplification steps are integrated in a single-tube workflow to generate the WGA product for library preparation.



**Figure 3. Takara Bio's patented library preparation method.** WGA products are enzymatically fragmented and ligated with stem-loop adapters, followed by UDI incorporation via PCR. Libraries are then ready for sequencing on Illumina platforms.

Following sequencing, cloud-based Embgenix Analysis Software (**Figure 4**) is used to analyze sequencing data from each sample. With Human genome build 38 (hg38) as an alignment reference, CNVs are identified and displayed in CNV plots as well as idiogram and list formats. Via a web-browser, users can easily perform CNV analysis, interactively review and amend CNV calls, and generate downloadable PDF reports with a variety of customizable features.


Input	Processing	Post-processing	Secondary-processing	Reporting
FASTQ files (results from Embgenix™ PGT-A Kit)	<ul style="list-style-type: none"> <li>Trim reads</li> <li>Align reads</li> <li>Count reads</li> </ul>	Compute calculated copy numbers (CCNs) using a proprietary algorithm	<ul style="list-style-type: none"> <li>Compute standard QC metrics</li> <li>Identify copy number variants (CNVs)</li> </ul>	<ul style="list-style-type: none"> <li>QC metrics</li> <li>CNV plot</li> <li>List of detected abnormalities</li> <li>Idiogram</li> </ul>

**Figure 4. Embgenix Analysis Software: general workflow for CNV detection and reporting**

## Results

### Detection of full chromosome and segmental aneuploidies in gDNA samples

The first study, using previously characterized gDNA samples, evaluated the feasibility and reliability of the Embgenix PGT-A Kit in detecting full-chromosome and segmental aneuploidies. The Embgenix software detected chromosomal abnormalities affecting whole chromosomes with high accuracy and confidence. An example plot resulting from this study is shown in **Figure 5**.



**Figure 5. Embgenix plot showing result for genomic DNA sample containing trisomy of Chr18.**

The Embgenix PGT-A Kit was applied to analyze 22 reference specimens with segmental aneuploidies. Embgenix PGT-A successfully detected 100% of the segmental aneuploidies, which ranged in size from 7.1 Mb to 91.4 Mb, with no false positives called. For 20 of the 22 segmental aneuploidies, Embgenix PGT-A provided highly accurate estimates (within 15%) of CNV size (**Table 1**). **Figure 6** shows example plots of accurate calls by Embgenix PGT-A for relevant CNV structures requiring high resolution of segmental detection (Coriell samples NA13019 and NA14485). In some rare cases, the CNV lies across a region of the genome where it is difficult to perform read mapping, the process of aligning reads to the reference genome. These “low mappability regions” may be excluded from the reported CNV size, resulting in significant underestimation of the reported CNV size. The underestimation of size for one CNV in Coriell sample NA20556 may have arisen due to this


issue, as that CNV extends into a low mappability region.

**Table 1. Study with 22 samples, showing accurate reporting of segmental aneuploidies including size prediction.**

Takara Bio's in-house verification			
Coriell Sample	Chr	True CNV Size (Mb)	Embgenix PGT-A Predicted Size (Mb)
NA20022	Chr3	61.1	64
NA14164	Chr13	41.2	48
NA12606	Chr13	41.2	42
NA09367	Chr6	35.2	35
NA11672	Chr10	26.2	26
NA10800	Chr4	25.2	26
NA05966	Chr14	21.2	21
NA11213	Chr2	17.5	18
NA09216	Chr2	16.7	17
NA08331	Chr13	12.1	13
NA10989	Chr9	12	12
NA14943	Chr2	7.8	9

Independent testing at Clinical Laboratory A			
Coriell Sample	Chr	True CNV Size (Mb)	Embgenix PGT-A Predicted Size (Mb)
NA13019	ChrX	91.4	98
NA13019	ChrX	56.8	59
NA14485	Chr8	31.1	34
NA06226	Chr16	21.3	22
NA05966	Chr14	21.2	22
NA06870	Chr18	15.4	15
NA20556	Chr15	12.3	7**
NA08331	Chr13	12.1	13
NA16362	Chr22	8*	8
NA14485	Chr8	7.1*	7

\*CNV size less than minimum specification for Embgenix Analysis Software. \*\*CNV occurred in region with low mappability.




**Figure 6. Analysis using the Embgenix PGT-A Kit demonstrates concordance with previously characterized segmental aneuploidy. Panel A.** Embgenix PGT-A correctly predicts two segmental aneuploidies in ChrX, linked to Turner Syndrome, in Coriell gDNA sample NA13019. **Panel B.** Embgenix PGT-A correctly predicts two segmental aneuploidies in Chr8, linked to hydrocephalus, in Coriell gDNA sample NA14485.

**Detection of mosaicism**

Next, we tested the accuracy of mosaic detection by Embgenix PGT-A. Two clinical laboratories, A and B, performed studies using mosaic samples created using mixtures of two samples with different aneuploidies. Laboratory A used mixtures of Coriell gDNA samples NA05966 (21.2 Mb trisomy in Chr14) and NA10925 (16.1 Mb monosomy in Chr7). The ratios tested are shown in **Figure 7**, which presents the results of the study as a comparison between Embgenix PGT-A and another NGS PGT-A method, Reproseq. Embgenix PGT-A successfully detected the

segmental mosaicisms represented in the gDNA mixtures at 50% and 75%. Reproseq did not detect the mosaicisms at 50% and 75%.



**Figure 7. Comparison of mosaic detection capability between Embgenix PGT-A and Reproseq.** Sensitivity of mosaic detection is compared side-by-side for Embgenix PGT-A and Reproseq, using artificial segmental mosaicism generated by mixing two Coriell gDNA samples (NA05966 and NA10925) containing segmental aneuploidies in different genomic locations. The two Coriell gDNA samples were mixed at the five ratios indicated. CNVs reported by each method are marked as true positives (blue arrows) or false negatives (red squares), with no false positives reported by either method.

The second Clinical Laboratory, B, performed a study using two gDNA samples with known CNVs. The first of these samples contained a 14 Mb segmental monosomy of Chr11. The second gDNA sample contained trisomy of Chr18 (78 Mb). Seven mixtures of these two gDNA samples were prepared. The resulting mixtures were processed using Embgenix PGT-A and Embgenix Analysis Software.

As shown in **Table 2**, Embgenix PGT-A successfully detected the Chr18 trisomy in all mosaic mixtures tested. Embgenix PGT-A successfully detected the 14 Mb Chr11 aneuploidy in mosaics where the sample with that aneuploidy made up at least 40% of the mixture.

**Table 2. Study of mosaic mixture demonstrating detection of full chromosome aneuploidy (78 Mb) and 14 Mb segmental aneuploidy.**

Embgenix PGT-A Kit validation using biopsy samples			
%mosaicism Chr11, 14 Mb segmental	%mosaicism Chr18, full chromosome	Embgenix PGT-A detection: Chr11, 14 Mb segmental	Embgenix PGT-A detection: Chr18, full chromosome
80%	20%	Yes	Yes
70%	30%	Yes	Yes
60%	40%	Yes	Yes
50%	50%	Yes	Yes
40%	60%	Yes	Yes
30%	70%	No	Yes
20%	80%	No	Yes

Multiple IVF clinics collaborated to validate the performance of Embgenix PGT-A with TE biopsy samples. The table in **Figure 8, Panel A**


summarizes two independent datasets generated in laboratories at two separate clinics, B and C. In the study performed by Clinical Laboratory B, 34 samples (including 26 TE samples and 8 gDNA samples) were analyzed using Embgenix PGT-A with Embgenix Analysis Software and VeriSeq PGS Kit with BlueFuse Multi Software. Embgenix results were 91% concordant with results obtained from the clinic's established NGS PGT-A analysis method. The three samples generating non-concordant results were from the same embryo. Clinical laboratory B determined these samples probably came from parts of the embryo with different genetic content.

In the study performed by Clinical Laboratory C, multiple biopsies from the trophectoderm (TE) or inner cell mass (ICM) of 71 embryos were analyzed using Embgenix PGT-A and other NGS PGT-A methods and software. Embgenix results were 100% concordant with results obtained from the established NGS PGT-A analysis for each embryo. **Figure 8, Panels B, C, and D** show a representative example comparing analysis of two TE samples and one ICM sample from Embryo 9. The results of analysis with Embgenix PGT-A were 100% concordant with data obtained by sample analysis with established methods and software. **Figure 9** shows results for another example from this study, Embryo 47.


**A**

Clinical Laboratory	Sample composition	Concordance
B	26 TE biopsies and 8 gDNA samples	91%
C	TE or ICM biopsies from 71 embryos	100%


**B**




**C**



**D**



**Figure 8. Aneuploidy assessments conducted by Clinical Laboratories B and C demonstrate concordance between results from Embgenix PGT-A analysis and results from established PGT-A methods.** Representative example data from Embryo 9 are shown in Panels B–C. Aneuploidy assessment of two TE biopsies and one ICM biopsy were performed using Embgenix PGT-A and an established method, VeriSeq PGS Kit with BlueFuse Multi Software. All analyses indicated monosomy of Chr5. **Panel A.** Table summarizing Embgenix PGT-A performance on datasets from collaborating IVF clinics. **Panel B.** Data from a TE sample assessed with VeriSeq PGS Kit with BlueFuse Multi Software. **Panel C.** Data from an ICM sample assessed with Embgenix PGT-A and Embgenix Analysis Software. **Panel D.** Data from a TE sample assessed with Embgenix PGT-A and Embgenix Analysis Software.



**Figure 9. Aneuploidy assessment conducted by Clinical Laboratory C demonstrates concordance between results from Embgenix PGT-A analysis and results from established PGT-A methods.** Representative example data from Embryo 47 is shown. Results from analysis of biopsy samples processed with Embgenix PGT-A and another PGT-A method are in concordance, indicating monosomy of Chr16. **Panel A.** Data from a TE sample processed with an established NGS PGT-A method and software. **Panel B.** Data from a TE sample processed with Embgenix PGT-A and Embgenix Analysis Software.

## Conclusion

We have shown that the Embgenix PGT-A Kit, by incorporating PicoPLEX WGA, Illumina library preparation, and Embgenix Analysis Software for copy number assessment, provides a complete solution for PGT-A.

Our study data demonstrate that Embgenix PGT-A rapidly and accurately calls whole chromosome and segmental aneuploidies greater than 10 Mb, as well as low (30–50%) and high (50–70%) mosaicisms for defined mixtures of samples with known aneuploidies.

The results also show that Embgenix PGT-A accurately determines the size of CNVs in segmental aneuploidies involving a wide array of different chromosomal regions. Testing in independent clinical laboratories using both gDNA and TE and ICM biopsies has shown concordance between results from Embgenix PGT-A and results from established approaches for NGS-based PGT-A.

Taken together, the results clearly show the feasibility of application of Embgenix PGT-A and establish the kit as a viable solution for current and future PGT-A applications.

## Materials and Methods

### Sample collection

For experiments performed by Takara Bio's in-house laboratory, cell lines with defined segmental aneuploidies were purchased from Coriell Biorepository (USA) as reported in **Table 3**. 26 TE biopsies were obtained by Clinical Laboratory B. A total of 71 embryos were included in the study conducted by Clinical Laboratory C. For each embryo, multiple biopsies were taken from the trophectoderm (TE) or inner cell mass (ICM).

### Table 3. Panel of Coriell gDNA samples containing segmental aneuploidies

Coriell ID	Description	Size (Mb)
NA05966	Duplicated chromosome	21.2
NA06226	Translocated chromosome	21.3
NA06870	Aneuploid chromosome number - non-trisomic	15.4
NA08331	Chromosome deletion	12.1
NA09216	Chromosome deletion	16.7
NA09367	Duplicated chromosome	35.2
NA10800	Chromosome deletion	25.2
NA10989	Gilles de la Tourette syndrome (GTS) Chromosome deletion	12.0
NA11213	Chromosome deletion	17.5
NA11672	Chromosome deletion	26.2
NA12606	Aneuploid chromosome number - non-trisomic	41.2
NA13019	Chromosome deletion Turner syndrome	91.4, 56.8
NA14164	Tetralogy of Fallot Chromosome deletion	41.2
NA14485	Inverted duplication deletion	7.1, 31.1
NA14943	Chromosome deletion	7.8
NA16362	Aneuploid chromosome number - trisomy	8.0
NA20022	Duplicated chromosome	61.1
NA20556	Isodicentric chromosome	12.3

### Whole genome amplification

Samples were processed following the Embgenix PGT-A protocol. The kit employs a single-tube protocol for the amplification of genomic DNA starting from samples stored in liquid N<sub>2</sub>, a standard practice for WGA starting material in this field.

### Sequencing

Sequencing-ready libraries were pooled and sequenced using a MiSeq® instrument with 24 samples per flow cell or using a NextSeq® instrument with 96 samples per flow cell.

### Genomic data analysis and aneuploidy screening

Sequencing data were analyzed using Embgenix Analysis Software.

## References

Kemper, J. M. et al. Preimplantation Genetic Testing for Aneuploidy: A Review. *Obstet Gynecol Survey* **74(12)**, 727–737 (2019).

Vendrell, X. et al. New protocol based on massive parallel sequencing for aneuploidy screening of preimplantation human embryos. *Systems Biology in Reproductive Medicine* **63(3)**, 162–178 (2017).



Zhang, X. et al. The comparison of the performance of four whole genome amplification kits on ion proton platform in copy number variation detection. *Bioscience Reports* **37** BSR20170252 (2017)

**Takara Bio USA, Inc.**

United States/Canada: +1.800.662.2566 • Asia Pacific: +1.650.919.7300 • Europe: +33.(0)1.3904.6880 • Japan: +81.(0)77.565.6999

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. © 2023 Takara Bio Inc. All Rights Reserved. All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at [takarabio.com](https://www.takarabio.com).