PrimeSTAR® GXL Premix Fast, Dye plus

Code No. R052S Size: 1 ml
(for 40 reactions)

Description:

PrimeSTAR GXL Premix Fast, Dye plus is an optimized 2X premix for fast PCR consisting of PrimeSTAR GXL DNA Polymerase (a high-fidelity DNA polymerase developed by Takara Bio), reaction buffer, dNTP mixture, gel loading dye (green), and a density reagent. PCR reactions performed with this product can be loaded directly into the gel for agarose gel electrophoresis. PrimeSTAR GXL Premix Fast, Dye plus is a robust, high-fidelity PCR enzyme that can be also used in difficult reactions, including amplification of long targets and GC-rich regions. This enzyme is designed for hot-start PCR using a monoclonal antibody that prevents non-specific amplification due to mis-priming and/or formation of primer dimers before the PCR reaction.

Storage: -20°C

Note: Freeze-thawing of reagents should be limited to ≤ 20 times. >20 freeze-thaw cycles may affect performance. Gently mix well and centrifuge briefly before use.

General Composition of PCR Reaction Mixture (total 50 μ l):

PrimeSTAR GXL Premix Fast, Dye plus (2X)	25 μ l (final conc. 1X)
primer 1	10 pmol (final conc. 0.2 μ M)
primer 2	10 pmol (final conc. 0.2 μ M)
Template	See (2) Template on the right
Sterile purified water	up to 50 μ l

PCR Conditions:

98°C	10 sec	} 30 cycles
55 or 60°C*1	5 sec	
68°C*2	1 - 10 sec/kb*3	

or

98°C	10 sec	} 30 cycles
68°C*2	1 - 10 sec/kb*3	

*1 When the Tm value of the primers is more than 55°C, set the annealing temperature to 60°C. When the Tm value is 55°C or less, set the annealing temperature to 55°C.

$$Tm (^{\circ}C) = [(the\ number\ of\ A\ and\ T) \times 2] + [(the\ number\ of\ G\ and\ C) \times 4] - 5$$

*2 For 3-step PCR, set the extension temperature to 68°C.

*3 Estimated extension time setting

amplification of products ≤ 10 kb:	5 sec/kb
amplification of products > 10 kb:	10 sec/kb

The yield of the amplification product may be improved by increasing the extension time.

Selecting PCR Conditions

- For amplification of products ≤ 10 kb in length, try to perform 3-step PCR at first.
- For amplification of products ≥ 10 kb in length, or of GC-rich targets 2-step PCR is recommended when enhanced specificity is desired.

Optimization of Reaction Parameters:

In order to obtain the best PCR results, it is important to optimize the PrimeSTAR GXL Premix Fast, Dye plus reaction parameters to fully utilize the enzyme's properties and advantages.

(1) Primer design

Select primer sequences using primer design software

[For ≤ 10 kb amplified products]

For general amplification, 20- to 25-mer primers are suitable. Selection of primers with a Tm value of $\geq 55^{\circ}C$ or greater than 25-mer in length may provide optimal results.

[For >10 kb amplified products]

Design primers that are 25- to 35-mers and that have a Tm value of $\geq 65^{\circ}C$. Avoid high GC-content at the 3' end of each primer.

[For GC-rich amplified products]

Design primers to have Tm values > 60°C.

Note: Do not use inosine-containing primers with PrimeSTAR GXL Premix Fast, Dye plus.

(2) Template

Recommended quantities of template DNA (for a 50 μ l reaction):

(for long PCR products)

Human genomic DNA	5 - 500 ng	(100 - 500 ng)
<i>E. coli</i> genomic DNA	100 pg - 200 ng	(10 - 200 ng)
Plasmid DNA	10 pg - 10 ng	(1 - 10 ng)
cDNA	25 - 750 ng	(250 - 750 ng)

Note: Do not use templates containing uracil, such as bisulfite-treated DNA.

Electrophoresis of Amplified Products:

TAE buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using PrimeSTAR GXL Premix Fast, Dye plus. The vivid green color separates into blue and yellow dye fronts when the PCR product is run on an agarose gel.

Note: TBE buffer may result in DNA band patterns that are enlarged at the bottom of the gel.

Termini of Amplified Products:

Most PCR products amplified with PrimeSTAR GXL Premix Fast, Dye plus have blunt-end termini. Therefore, they can be cloned directly into blunt-end vectors. If necessary, phosphorylate amplified products before cloning. Use the Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) for cloning into a blunt-end vector.

If you want to clone amplified products into T-vectors, dA addition to the 3' end is necessary.

Restriction Enzyme Digestion:

Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of PrimeSTAR GXL Premix Fast, Dye plus from the reaction mixture by phenol/chloroform extraction or by using the NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/.50/.250)* kit. A 3'-overhang produced with restriction enzyme, such as *Pst* I, may be removed by the 3' \rightarrow 5' exonuclease activity of PrimeSTAR GXL Premix Fast, Dye plus if residual polymerase remains in the restriction enzyme reaction.

* Not available in all geographic locations. Check for availability in your area.

PrimeSTAR is a registered trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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