# PrimeSTAR® LongSeq DNA Polymerase

**Code No. R055S**        **Size:            1 ml  
  (for 40 reactions)**

**Code No. R055A**        **Size:            1 ml x 5  
  (for 200 reactions)**

## Description:

PrimeSTAR LongSeq DNA Polymerase is a novel PCR enzyme based on DNA polymerase derived from *Thermococcus* sp. This product improves the PCR success rate for GC or AT rich sequences and long targets that are difficult to amplify with general PCR reagents due to the optimization of enzymes and buffers. Furthermore, its superior specificity allows for long-range multiplex PCR. This product is 2X premix type that contains all components (enzyme, buffer, dNTP Mixture, etc.) and does not freeze at -20°C, so it is possible to perform PCR by adding primers and templates without thawing. This enzyme is designed for hot start PCR using the monoclonal antibody. Thus, it is possible to prevent non-specific amplification due to mispriming and/or formation of primer dimers before the PCR reaction.

## <Features>

- High fidelity
- Exceptional specificity
- Capable of long-range multiplex PCR, which was challenging with conventional products
- Convenient 2X premix type

**Storage:** -20°C

Note: Gently mix well before use and centrifuge briefly.

## Electrophoresis, Cloning, and Sequencing of Amplified Products:

- (1) Electrophoresis of amplified products:  
TAE Buffer is recommended for agarose gel electrophoreses of PCR products amplified using PrimeSTAR LongSeq DNA Polymerase. TBE Buffer may distort bands, particularly small products at the bottom of the gel.
- (2) Termini of amplified products:  
Most PCR products amplified with PrimeSTAR LongSeq DNA Polymerase have blunt-end termini. Accordingly, they can be cloned directly into blunt-end vectors.  
PCR products amplified with primers designed for in-Fusion cloning can be cloned as is. If necessary, phosphorylate amplified products before cloning. In-Fusion® Snap Assembly Master Mix (Cat. #638947/638948/638949/638943/638944/638954/638955/638956) can be used for In-Fusion cloning. Use the Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) for cloning into a blunt-end vector.  
If you want to clone amplified products into T-vectors, dA addition to the 3' end is needed.
- (3) Restriction enzyme digestion:  
Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of PrimeSTAR LongSeq DNA Polymerase from the reaction mixture by phenol/chloroform extraction or by using NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/.50/.250)\*. Particularly for 3'-protruding restriction enzymes such as *Pst* I, the 3'-protruding termini produced by these enzymes may be deleted by 3' → 5' exonuclease activity of PrimeSTAR LongSeq DNA Polymerase if residual polymerase remains present in the restriction digest reaction.
- (4) Direct sequencing  
Perform phenol/chloroform extraction of PCR products prior to direct sequencing to ensure inactivation of 3' → 5' exonuclease activity.  
Alternatively, NucleoSpin Gel and PCR Clean-up (Cat. #740609.10/.50/.250)\* or NucleoSpin gDNA Clean-up (Cat.#740230.10/.50/.250)\* may be used to purify DNA prior to sequencing.

\* Not available in all geographic locations. Check for availability in your area.

## General Composition of PCR Reaction Mixture (Total 50 µl):

PrimeSTAR LongSeq Premix (2X)	25 µl (final conc. 1X)
primer 1	10 pmol (final conc. 0.2 µM <sup>*1</sup> )
primer 2	10 pmol (final conc. 0.2 µM <sup>*1</sup> )
Template	See "Template" chapter
Sterile purified water	up to 50 µl

## PCR Conditions:

94°C	1 min*2	
↓		
98°C	10 - 20 sec	] 30 - 40 cycles
55 - 67°C*3	15 sec	
68°C	20 - 30 sec/kb	
or		
98°C	10 - 20 sec	] 30 - 40 cycles
68°C	20 - 30 sec/kb	

- \*1 For multiplex PCR, it is recommended to prepare each primer to a final concentration of 0.2 µM.
- \*2 When amplifying a GC-rich region or a long target, perform initial denaturation at 94°C for 1 min.
- \*3 Please set the annealing temperatures according to the  $T_m$  value of the primers.

Note 1: Do not use inosine-containing primers with PrimeSTAR LongSeq DNA Polymerase.

Note 2: Ensure that template and primers do not contain excessive EDTA.

Note 3: TE buffer (10 mM Tris-HCl, 0.1 mM EDTA) or sterile purified water is recommended for the dilution of template and primers.

## <Primer Design>

- Design primers that are approximately 35-mers with a  $T_m$  value of > 70°C.
- If the  $T_m$  value is below 70°C, consider using 3-step PCR for better amplification.
- Especially for multiplex PCR, it is recommended to design primers so that the  $T_m$  values of each primer are approximately equivalent.

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PrimeSTAR is a registered trademark of Takara Bio Inc.  
In-Fusion is a registered trademark of Takara Bio USA, Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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