BcaBEST™ DNA Polymerase ver.2.0

Code No. RR380A **Size:** **1,600 U**
Conc.: **8 U/μl**

Supplied Reagents:

BcaBEST DNA Polymerase ver.2.0 **200 μl**
2X BcaBEST Buffer **1.25 ml x 2**

Description:

BcaBEST DNA Polymerase ver.2.0 is an updated product of the BcaBEST DNA Polymerase that lacks 5' → 3' exonuclease activity, originating from a thermophilic bacterium *Bacillus caldotenax*. By introducing mutations into the protein, the activity of both DNA polymerase and reverse transcriptase is highly enhanced while exerting the original strong strand-displacement activity. Therefore, the isothermal amplifications using this enzyme enable to amplify the target product from a little nucleic acid in a shorter time. This product also contains a reaction buffer including dNTP, which is directly used for LAMP (Loop-Mediated Isothermal Amplification) assays for testing the detection of a certain target gene.

Storage: -20°C

Source:

Escherichia coli carrying a plasmid containing the gene for BcaBEST DNA Polymerase ver.2.0

Properties:

Molecular mass: approx. 68 kDa

Unit definition:

One unit is defined as the amount of enzyme that incorporates 210 nmol of dNTP into acid-insoluble precipitate in 30 minute at 65°C.

Reaction Mixture for Unit Definition:

25 mM	TAPS, pH9.3 at 25°C
50 mM	KCl
2 mM	MgCl ₂
0.1 mM	DTT
200 μM	dATP • dGTP • dCTP each
100 μM	dTTP
10 μCi	[³ H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. LAMP (Loop-Mediated Isothermal Amplification)
2. RCA (Rolling Circle Amplification)

Precautions for use:

1. Don't mix the enzyme vigorously.
2. Extreme cautions should be taken during reaction preparation to avoid contamination. If analyzing the products after reaction, execute in a different area from the reaction preparation.
3. Reaction temperatures above 70°C are not recommended. Hence this enzyme can't be used for conventional PCR.

LAMP Reaction Example:

Nuclease-free water	x μl
2X BcaBEST Buffer	12.5 μl
10X LAMP primer mix*1	2.5 μl
BcaBEST DNA Polymerase ver.2.0 (8 U/μl)	0.3 - 1 μl*2
Sample*3	e.g. 2 μl
<hr/>	
Total	25 μl

Incubate at 60 - 65°C for 30 minutes.*4

- *1 FIP/BIP primers : 16 μM
 F3/B3 primers : 2 μM
 LoopF/B primers : 8 μM

*1,2,4 Optimize reaction systems, including primer and enzyme concentrations, and reaction temperature and time, due to their variations depending on the target and primer sequence. Use of high concentration of enzyme would reduce the detection (amplification) time but might also increase the risk of non-specific amplification. It is strongly recommended to include a no-template control reaction to ensure the specificity.

*3 In order to prevent the non-specific reaction during its preparation, sample should be added lastly. In case of RNA target, add the Reverse Transcriptase XL (AMV) (5 U/μl) to be a final concentration of 0.04 U/μl.

Note: Various detection methods can be combined with the reaction, such as a color indicator detection by eye and a fluorometric detection by real-time instrument. Use an appropriate method depending on your purpose.

References:

- 1) Tomita N, Mori Y, Kanda H, and Notomi T. *Nat Protoc.* (2008) **3**: 877-882.
- 2) Zhao Y, Chen F, Li Q, Wang L, and Fan C. *Chem Rev.* (2015) **115**: 12491-12545.

Related products:

Reverse Transcriptase XL (AMV) for RT-PCR Kit (Cat. #2630A)
RNase-free Water (Cat. #9012)

BcaBEST is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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